




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



Technical sheet

Table Cover

Reference

1900CEB, 2100CEB, 8186CEB, 8368CEA, 8371CEB, 8378CEB, 8380CEB

Description

Reference	Size (cm)	Reinforcement width (cm)	Thickness (µm)	Unit per case	Sterile wrap dimensions (in cm)	Packaging
1900CEB	140x190	74	75	22	60x60	Peel pouch
2100CEB	140x120	66	75	26	60x60	Peel pouch
8186CEB	160x240	74	90	14	60x60	Header bag
8368CEA	120x140	66	50	30	60x60	Peel pouch
8371CEB	140x190	74	50	22	60x60	Peel pouch
8378CEB	125x230	74	50	18	60x60	Peel pouch
8380CEB	203x240	74	75	14	76x76	Header bag

Main materials

Polyethylene
SPP reinforcement: Absorbent Polypropylene (spunbond)

Packaging

Peel pouch or Header bag
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795 standard



© 2012 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-1900CEB-2100CEB-8186CEB-8368CEA
8371CEB-8378CEB-8380CEB-EN03

01.01.INT.12