



Product description

Order number	95.64.997
Product description	Storage box made of cold-resistant, laminated cardboard, with lid and compartmentalized insert, for tubes up to max. \varnothing 12 mm and a tube length of 36-45 mm, grid: 10 x 10, for 100 tubes, External dimensions (WxLxH): 135 x 135 x 45 mm, 1 pc./bag

Product characteristics

Cap	Slip-on lid
Suitable for	Tubes up to max. 12 mm diameter, tube height between 36-45 mm
Storage capacity	100
Format	10 x 10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



Size

Width of product	135 mm
Length of product	135 mm

Material & colours

Product material	Cardboard
Cap material	Cardboard
Colour of cap	White

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	batched

Packaging

Minimum order qty.	1
Pcs. / inner box	1
Pcs. / box	120
Pcs. / pallet	1440
Case volume	0,1139 cbm
Weight of case	11,5 KG
EAN of inner box	4038917095907
EAN of outer box	4038917007887

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

