







SPIRO

	REF.	DESCRIPTION	QTY
	121800	<p>Spiro phonation valve set consisting of:</p> <ul style="list-style-type: none"> • Housing with phonation valve and filter • Housing with filter, without phonation valve • Deflector • Connectors for silver cannulas • Tweezers <p>• With oxygen port</p> <p>• Single patient use</p>	1
	121801	<p>Spiro filter housing with phonation valve</p> <ul style="list-style-type: none"> • Single patient use 	5
	121802	<p>Spiro filter housing without phonation valve</p> <ul style="list-style-type: none"> • Single patient use 	5
	121803	<p>Spiro filter for Spiro filter housing 121800/121801/121802</p> <ul style="list-style-type: none"> • Single use 	100

Soft Membrane Phonation Valve

	REF.	DESCRIPTION	QTY
	121450	<ul style="list-style-type: none"> • For 15 mm connectors • Soft Membrane technology and O₂-adapter • Single patient use 	2

Cleaning Brushes

	REF.	DESCRIPTION	QTY
	968675	<ul style="list-style-type: none"> • Size 6 mm for tube I.D. 3–6 mm • Size 10 mm for tube I.D. 7–9 mm • Size 15 mm for tube I.D. 10–13 mm 	10