



PRODUCT SPECIFICATION

Lid 69 white



PRODUCT

Product name: Lid 69 white
Article number: C112040
Description: Fits in a container 170 ml and other

DESIGN

Dimensions:
External height: 6.2 mm
Outer diameter: 70.5 mm
Weight: 3.9 g

MATERIAL

Raw material: Polyethylene CA 9150
Masterbatch: Masterb vit LDPE 1001096

PRODUCTION

Process: Injection moulding.
Hygiene: Production and packing in hygienic environment

PACKAGING

2.200 lids are packed into a LDPE bag.
The LDPE bag is placed in a cardboard outer carton
(570 x 380 x 410).

Labeling: Each box is labeled with the following information:

Product name:
Production date:
Article number:
Quantity:
Batch number