


Bone & Bone Marrow Biopsy

T-Lok™ Bone Marrow Biopsy Needles

T-Lok Bone Marrow Biopsy Needles are designed to provide quality biopsies by entrapping bone marrow in the needle without altering the sample.



Box of 10

Catalogue #	Gauge	Length	Cannula
DBMNJ0804TL	8G	4" (10.2cm)	w/Extraction Cannula
DBMNJ0806TL	8G	6" (15.2cm)	w/Extraction Cannula
DBMNJ1104TL	11G	4" (10.2cm)	w/Extraction Cannula
DBMNJ1106TL	11G	6" (15.2cm)	w/Extraction Cannula
DBMNJ1104D	11G	4" (10.2cm)	w/D Shape Extraction Cannula