

RHINOLOGI

METALLSUG MEDIPLAST® TYP NASAL

6066516020	2,0 x 150 mm	grön	luer-lock
6066516030	3,0 x 150 mm	blå	rund hub
6066516040	4,0 x 150 mm	gul	rund hub
6066502150	2,0 x 150 mm	grön	rund hub

METALLSUG MEDIPLAST®

Med transparent handtag och sugkontroll

6066509120	2,0 x 110 mm	6066516120	2,0 x 150 mm
6066509125	2,5 x 110 mm	6066516130	3,0 x 150 mm
6066509130	3,0 x 110 mm	6066516140	4,0 x 150 mm
6066509135	3,5 x 110 mm		

Med transparent handtag, 2 sidohål och sugkontroll

6066509020	2,0 x 110 mm
6066509025	2,5 x 110 mm
6066509030	3,0 x 110 mm
6066509035	3,5 x 110 mm

METALLSUG MEDIPLAST® TYP FERG

Med transparent handtag och sugkontroll

6066501020	2,0 x 150 mm
6066501030	3,0 x 150 mm
6066501040	4,0 x 150 mm

Som ovan men med sidohål

6066511020	2,0 x 150 mm
6066511030	3,0 x 150 mm
6066511040	4,0 x 150 mm

METALLSUG MEDIPLAST® TYP FRAZ

Med transparent handtag och sugkontroll

6066508020	2,0 x 110 mm
6066508030	3,0 x 110 mm
6066508040	4,0 x 110 mm

METALLSUG MEDIPLAST® TYP ENDO


Med transparent handtag utan sugkontroll

6066503020	2,0 x 150 mm
6066503025	2,5 x 150 mm
6066503030	3,0 x 150 mm
6066503035	3,5 x 150 mm
6066503040	4,0 x 150 mm

METALLSUG MEDIPLAST® TYP ZOELLNER

Med transparent handtag

6066504020	2,0 x 130 mm	utan sugkontroll
6066504025	2,5 x 130 mm	utan sugkontroll
6066504030	3,0 x 130 mm	utan sugkontroll
6066514025	2,5 x 90 mm	med sugkontroll I



typ Nasal



transp handtag



transp handtag



typ Ferg



typ Fraz



typ Endo



typ Zoellner