



# Extension set

MF2299E



- Needle Free System
- Set Length 205 cm
- Priming Volume 1 ml
- Bore 0.7 mm
- Flow Rate 0 - 750 ml/hr
- 1 Pinch Clamp(s)
- Pressure Sensor Disc
- 1 Male Luer(s)
- Rotating Male Luer(s)
- Low Sorbing
- Polyethylene lined PVC
- Latex-Free
- Quantity 100



Customer Service  
0800 917 8776 freephone  
+ 44 (0)1256 388 284  
(0)1256 330 860

CareFusion UK 306 Ltd  
Basingstoke