

PRODUKTINFORMATION TRAKEOSTOMIKANYL



UNIPERC®, PVC

- För vuxna med avvikande anatomi eller vävnadsvullnad eller "XL-patienter"
- Justerbar kanyl
- Mjuk och armerad termosensitiv PVC
- Stor mjuk krage med information om ID och OD
- Soft Seal® högvolum lågtryckskuff*
- Pilotballongen är försedd med information om storlek och kuffens vilodiameter
- 15 mm ISO standard konnektor
- Tydlig längdmarkering på kanylen
- Mjuka avfasade kanter vid kanylens ände
- Finns med och utan Soft Seal® kuff
- Tunn flexibel innerkanyl som inte minskar trakealkanylens innerlumen på grund utav specialdesign



*Pilotballongen ska placeras och fixeras utanför det område som ska scannas.

Benämning	Art nr	Enhet	Avdfp
Uniperc® justerbar trakeostomikanyl med kuff stl 7,0 mm	100/897/070	st	1
Uniperc® justerbar trakeostomikanyl med kuff stl 8,0 mm	100/897/080	st	1
Uniperc® justerbar trakeostomikanyl med kuff stl 9,0 mm	100/897/090	st	1
Uniperc® justerbar trakeostomikanyl utan kuff stl 7,0 mm	100/899/070	st	1
Uniperc® justerbar trakeostomikanyl utan kuff stl 8,0 mm	100/899/080	st	1
Uniperc® justerbar trakeostomikanyl utan kuff stl 9,0 mm	100/899/090	st	1

Extra tillbehör;

- Extra guidewire och innerkanyler kan beställas separat
- Kan levereras med speciell designat Uniperc® PDT set

Smiths Medical Sverige AB

Box 1143, 164 22 Kista
Tlf. 08 - 594 772 50
Fax.08 - 594 772 69
info.sweden@smiths-medical.com
www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
©2012 Smiths Medical. All rights reserved. SE_13cGPI_100897+9_1