




1.5 litre USON Sharps Container

This USON 1,5 litre Sharp Containers is designed for hospitals and health care institutions to dispose of syringes, scalpels and other types of sharps.

It is designed to maximise safety whilst still maintaining a barrier to protect the user and prevent spillage during use. Different lid designs give the possibility to discharge needles, syringes and other types of sharps.

USON sharps containers meet the UN requirement no. 3291 regarding packaging approved for transportation of dangerous goods (waste/disposal)

- Products are made of environmentally friendly, strong and durable Polypropylene
- Design and manufacturing secure a high degree of safety
- Different designs of lid are provided for different applicability
- Reversible closure can be safely fixed
- Containers are waterproof when closed
- Standard colours are yellow and yellow transparent
- Products are made with a private label and own colour when ordered in large quantities
- Containers can be supplied without UN certificate (for less costly solutions)

MATERIAL

Polypropylene (PP)



DIMENSIONS

Diameter (top): 13,70 cm

Height: 13 cm

PACKAGING

Per carton: 100 pcs.

Per pallet: 2.000 pcs.

ACCESSORIES

Disc to fix on table or trolley

Hangers for mounting on wall and trolley

Disc #4799



Wall hanger #5200 Trolley hanger #5200 + #6000

ISO 9001 · ISO 13485
Management System Certification
BUREAU VERITAS
Certification Denmark A/S



USON PLAST A/S is
certified by BVQI
according to DS/EN
ISO 9001 and 13485