

PRODUKTINFORMATION

PRO-VENT® PLUS MED EDGE™ & FILTER-PRO®



BLODGASSPRUTA MED KANYL

Blodgassprutan Pro-Vent® Plus med sticksäker kanyl, EDGE™ och Filter-Pro® är konstruerad för provtagning av blodgaser. Pro-Vent® Plus har ett självluftande membran. Frystorkat elektrolytbalanserat litiumheparin är sprayat på insidan av sprutcyklindern vilket underlättar blandning.

- Självfyllande till önskad förinställd volym, går även att aspirera.
- EDGE™ nålskydd tillåter arbete med enhandsfattning och bibehållen provtagningsteknik.
- En pilindikator visar nålsögats position.
- Filter-Pro® är en skyddshylsa som ger en enkel och säker eliminering av luft utan risk för blodkontakt.
- Filter-Pro® fungerar även som skyddhatt vid transport av blodgassprutan.



Benämning	Färgkod	Nål-storlek (mm)	Sprut-storlek (mm)	Filter-Pro®	Luer lock	Luer slip	Art nr	Enhet	Avdfp
EDGE™ Pro-Vent® Plus	Orange	25G x 16	1	Ja	Nej	Ja	4619PG	st	200
EDGE™ Pro-Vent® Plus	Blå	23G x 25	1	Ja	Nej	Ja	4629PG	st	200
EDGE™ Pro-Vent® Plus	Grå	22G x 25	3	Ja	Nej	Ja	4645PG	st	200
EDGE™ Pro-Vent® Plus	Blå	23G x 25	3	Ja	Nej	Ja	4698PG	st	200
EDGE™ Pro-Vent® Plus	Blå	23G x 25	3	Ja	Ja	Nej	4558PG	st	200
EDGE™ Pro-Vent® Plus	Grå	22G x 25	3	Ja	Ja	Nej	4590PG	st	200

Material; Polypropylene, Acrylic och polyetylen filter, medicinsk silikonolja, rostfritt stål, frystorkat elektrolytbalanserat litiumheparin

Smiths Medical Sverige AB

Box 1143, 164 22 Kista
 Tlf. 08 - 594 772 50
 Fax.08 - 594 772 69
 info.sweden@smiths-medical.com
www.smiths-medical.com

Smiths Medical is part of the global technology business Smiths Group plc. Product(s)described may not be licensed or available for sale in all countries. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
 ©2012 Smiths Medical. All rights reserved.SE_11oGPI_PV_EDGE-2