




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



ECLIPSE®

Hand drape

Reference(s)

29419CEA

Description

A sterile wrap 86x86 cm contains:

- 1 spunlace hand drape 196x371 cm, 45x86 cm reinforcement, armboard cover, 3 cm diameter elastic fenestration, 4 cord holding tabs

Main materials

Eclipse® spunlace (polyester + woodpulp)

Polyrayon absorbent reinforcement (polyester + viscose)

Pack factor

10 units per carton

Packaging

Vented sterilization pouch

1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS

EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC

In compliance with the standard EN13795

MDD Class: I sterile



Latex free



0086



STERILE EO

© 2013 Medline Industries Inc. or one of its subsidiaries.
All rights reserved..

LIT-29419CEA-EN03

01.01.INT.12