



## Argon Safedraw™-P Safe Arterial Sampling DTX Plus™ Disposable Transducer Sets

No. of disposable transducers	Flush/Flow Rate (cc/Hr)	Total Pressure Line Length (cm)	Pressure Line Colour	Additional Information	Safedraw Syringe Volume (cc)	Septum Model	Box/Case Content (Units)	Product Code
1	3	183cm (153/30)	Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw™-P Septum (Red)	2cc	1 x Inline	5/20	688413
1	3	213cm (183/30)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw™-P Septum (Red)	2cc	1 x Inline	5/25	688410
1	3	223cm (133/90)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw™-P Septum (Red)	3cc	1 x Inline	5/20	688457
1	3	78cm (63/15)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw™-P Septum (Red) with Velcro Strap for Arm mounting	1.5cc	1 x Inline	5/20	688538
1	3	A: 213cm (183/30) V: 150	A: Red/Clear V: Blue	Single Admin Set. Single Transducer Set with non-shunting stopcock and Integrated Safedraw™-P Septum (Red) (with Separately Packaged 150cm Blue pressure line)	2.5cc	1 x Inline	5/25	688403