

Product Identification

Product family	Sterile suture sets
Reference	10608032
Name	Suturset med metallinstrument
Sterilization	Ethylene oxide
Brand	Sterisets®



Material

Code	Designation	Qty	Composition
1006	Wrapping field PE/Tissue 50 x 70 cm green	1	Polyethylene/tissue 10 µm
219023INT	Kidney bowl PP 300 ml 214 x 120 x 43 mm blue w. grad.	1	Polystyrene
217018INT	Bowl PP 250 ml Ø 101 mm height 54 mm blue w. grad	1	Polypropylene
56502	Ball NW 13 x 15 cm Ø 35 mm 30 g/m2 with ring	10	Non woven 30 g/m2
570016	Scissors suture S/B STR 14,5 cm w/c	1	AISI-420 Stainless steel
570034	Forceps dressing delicate Adson 1 x 2 teeth smooth 12 cm w/c	1	AISI-410 Stainless steel
570021	Needle holder Mayo-Hegar X-SERR jaw 15 cm w/c	1	AISI-410 Stainless steel

Packaging

Type	Sterile single pack
Package weight	170g
Package dimensions	300 x 200 x 60mm
Quantity per inner box	NA
Quantity transport carton	42pcs

Hazards Identification

Latex	No
Medicinal substances	No
Animal tissue	No
Toxic or dangerous materials	No

Shelf life

5 years from the date of manufacture, unless any component has a shorter shelf-life.

Storage conditions

Keep away from sunlight and keep dry.

Environmental and disposal considerations

Not regarded as dangerous for the environmental. Dispose of material according to local, regional and national regulations.

Compliance Regulations & Standards

Regulation (EU) 2017/745 (as applicable), Directive 93/42/EEC, ISO 13485, ISO 14971, ISO 15223-1, ISO 10993-1, ISO 11737.

