

Fisherbrand™ Traceable™ Nano™ Timer

Offers ultra-small size with intuitive operation

Brand: Fisherbrand™



Description

- Counts down or up to 23 hours, 59 minutes, 59 seconds on two timing channels
- Accuracy: 0.01%
- 0.375 in. high digits may be read from 10 ft. (3m)
- For repetitive tests, two-channel memory returns display to original programmed time

Order Info

Dimensions: 2.5 x 2 x 0.5 in. (6.4 x 5.1 x 1.3cm); Weight: 2.5 oz. (71g)

Certifications

Individually serial-numbered Traceable Certificate is provided from an ISO 17025 calibration laboratory accredited by A2LA™. It indicates traceability to standards provided by NIST™

Specifications

Accuracy	0.01%	Alarms	High db alarm sounds for 1 min.
Certifications/Compliance	ISO 17025, A2LA, NIST	Dimensions (L x W x H)	2.5 x 0.5 x 2 in. (6.4 x 1.3 x 5.08cm)
Height (Metric) Display	0.93cm	Height (Metric)	5.08cm
Includes	Traceable certificate	Length (Metric)	6.4cm
Resolution	1 sec.	Time Range	23 Hours, 59 min., 59 sec.
Weight (Metric)	71g	Width (Metric)	1.3cm
Display	LCD	Product Type	Timer