


Förlängningsslang MaxPlus® Clear

MP9001C-0006



- Nålfritt system
- 20 cm Aggregatlängd
- Fyllnadsvolym 0.5 ml
- Bore 1.2 mm Diameter
- 1 Skjutklämma
- Roterande han-luer
- Latexfri
- Kvantitet 50



50



Latexfri



DEHP-fri



86

Customer Service
+ 46 (0) 8 5444 3200
(0) 8 5444 3225
CS.Sweden@carefusion.com

Carefusion Sweden 314 AB
Sollentuna