



# Alaris™ VP Series DEHP FREE Low Sorbing Infusion Set with 15µm chamber filter

70953



- 300 cm Aggregatlängd
- Fyllnadsvolym 27 ml
- Bore 3 mm Diameter
- 1 Droppkammare med spike (ventilerad(e)/oventilerad(e))
- Hydrofobiskt luftfilter (med lock)
- Droppar per ml 20
- 1 Rullklämma
- Pumpsegment
- Alaris® säkerhetsklämma
- Tryckavkännare
- 1 Han-luer
- Lågabsorberande
- Latexfri
- Kvantitet 80



Customer Service  
+ 46 (0) 8 5444 3200  
(0) 8 5444 3225  
[CS.Sweden@carefusion.com](mailto:CS.Sweden@carefusion.com)

Carefusion Sweden 314 AB  
Sollentuna