



Enter Search Item

- Shop
- Support
- Repair
- Learn

Sign In | Register

Home > Shop > Accessories > Anesthesia Delivery > Gas Monitoring > D-FEND WATER TRAP, BLACK, PKG OF 10 PCS

## D-FEND WATER TRAP, BLACK, PKG OF 10 PCS



Item No.	876446-HEL
Product Category	Anesthesia Delivery
Product Type	Gas Monitoring

[Add to Wish List](#)

Login To Shop

[Product Specific Terms](#)



### Description:

### Compatible Products

D-FEND WATER TRAP, BLACK, PKG OF 10 PCS

SHOP	SUPPORT	REPAIR	LEARN	MY ACCOUNT
<ul style="list-style-type: none"> <li>Parts</li> <li>Accessories</li> <li>Equipment</li> </ul>	<ul style="list-style-type: none"> <li>Service Documents</li> <li>Lifecycle Information</li> <li>Part Identification Tool</li> </ul>	<ul style="list-style-type: none"> <li>Portable Device Repair</li> <li>Ultrasound Probe Repair</li> <li>Track My Repairs</li> </ul>	<ul style="list-style-type: none"> <li>How to Videos</li> </ul>	<ul style="list-style-type: none"> <li>Account Profile</li> <li>My Order History</li> <li>Manage Members</li> </ul>

GE Healthcare | Privacy Statement | Terms and Conditions | Accessibility | Cookie Notice & Disclosure