



Product Specification

**Multi-purpose container, max. volume: 500 ml,
(LxØ): 150 x 70 mm, graduated, PP, transparent**

Page 1



Product description

Order number	75.9922.813
Product description	Multi-purpose container, max. volume: 500 ml, (LxØ): 150 x 70 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, sterile, 40 piece(s)/bag

Product characteristics

Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	500 ml

Size

Sample volume	500 ml
Max. volume	500 ml
Diameter	70 mm
Height of product	150 mm
Length of product	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-09-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Multi-purpose container, max. volume: 500 ml,
(LxØ): 150 x 70 mm, graduated, PP, transparent**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

Purity & certification

Product category	not regulated
CE certified	no CE mark
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	80
Type of smallest subpackaging	bag
Piece(s) / inner box	40
Piece (s) / outer case	80
Piece(s) / pallet	1280
Depth of case	548 mm
Width of case	373 mm
Height of case	407 mm
Case volume	0.0832 cbm
Weight of product	0.0655 kg
Weight of case	5.27 kg
EAN of inner box	4066916011338
EAN case	4066916011345

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-09-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT