

CHIRANA® 3-part disposable syringes



According to EN, ISO standards | Quality management certified according to EN ISO 13485:2016 | CE marked

- ✓ ergonomic design, smooth flow, highly transparent barrel, tight
- ✓ black extended graduation, easy to read
- ✓ designed tip for best connection, scratch-free for glove manipulation
- ✓ available in Luer and Luer lock
- ✓ made of polypropylene and lubricated synthetic rubber
- ✓ pyrogen-free, latex-free, DEHP-free, PVC-free, non-toxic
- ✓ sterile, individually blister packed
- ✓ no dead space volume for 1ml
- ✓ centrally positioned Luer lock

CHIRANA® 3-part disposable syringes

Size	Pieces per Box	Pieces per Case	Item Code LUER	Item Code LUER LOCK
1 ml	100	2000	CH03001L	CH03001LL
2 ml	100	3000	CH03002L	CH03002LL
3 ml	100	3000	CH03003L	CH03003LL
5 ml	100	2400	CH03005L	CH03005LL
10 ml	100	1600	CH03010L	CH03010LL
20 ml	50	800	CH03020L	CH03020LL
30 ml	50	400	CH03030L	CH03030LL
50 ml	---	100	CH03050L	CH03050LL