

Alarm thermometer, MM125

Leverantör: DICKSON




MM125 is a simple, cost effective temperature monitoring tool which is used as a backup to a traditional data logger in vaccine or drug monitoring application.

- ▶ Tamper-resistant audible/visual alarms
- ▶ Two thermistor probes submersed in glass beads for stabilized readings
- ▶ Tracks min/max temperatures
- ▶ Easy to read display

Leveransinformation: The MM125 supplied with the temperature sensor submersed in glass beads instead of a glycol solution (glass beads and glycol are recommended by the CDC for vaccine (including VFC) storage.)

ORDER

VWR art.nr	Antal	Benämning	Mätområde	Noggrannhet	Upplösning	B×D×H	Vikt	Givare	För
DICKMM125		Digital alarm thermometer for vaccine and drug monitoring	-50...+50 °C	±1,0 (-18...+49) °C	0,1 °C	46×27,5×8 mm	90,72 g	2 thermistors	1