



Återförsäljare:

**MEDI PLAST**

Mediplast AB

Box 1004, 212 10 Malmö

T 020-78 80 35

mediplast.info@mediplast.com

www.mediplast.com



## Technical Sheet

# Absorbent cotton OR towels

### Reference(s)

References	Qty/Pack	X-Ray detection	Loop	Colour	Unit/Case
28200004CEA	4	No	No	Green	20
28700002CEA	2	No	No	Blue	40
28700004CEA	4	No	No	Blue	20
28700002XR	2	Yes	Yes	Blue	40
28700004XR	4	Yes	Yes	Blue	20

### Description

References	Dimensions	Weight of towel	Number of layers when folded
28200004CEA	66x43 cm	60g+/-5g/piece	12 layers
28700002CEA	66x43 cm	60g+/-5g/piece	12 layers
28700004CEA	66x43 cm	60g+/-5g/piece	12 layers
28700002XR	66x43 cm	60g+/-5g/piece	12 layers
28700004XR	66x43 cm	60g+/-5g/piece	12 layers

### Main materials

- Towels are made of 90 % cotton fabric + 10 % viscose
- X-ray made of Barium sulfate, the Barium content  $\geq 60\%$
- The loop is stitched

### Packaging

Towels wrapped in paper. Paper wrapped towels packed in a pouch.

### Manufacturer

Medline International France SAS

EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC

MDD Class: I Sterile



STERILE EO

© 2017 Medline Industries Inc. or one of its subsidiaries. All rights reserved.

TDS\_COTTON\_OR\_TOWELS\_EN06

01.01.INT.12