

microlife®

A partner for people. For life.



Atrial Fibrillation detection with

AFIBsens

The first step in stroke prevention.



Product video



Bluetooth®

BP B6 Connect Blood Pressure Monitor

- PC connection with BP Analyzer software and USB cable included
- Comfortable, preformed cuff for most arm sizes (22–42 cm)
- MyBP: diagnostic 7-day average following international hypertension guidelines

Key features



High accuracy in Atrial Fibrillation (AF) detection for stroke prevention



IHB for irregular heartbeat detection



Gentle+ ensures optimal speed and pressure control for comfortable measurements



LED traffic light indicator for reading blood pressure classification



Simple arrow display to compare measured blood pressure with stored average



2 users with 99 data memory per user



Cuff fit check



Microlife Connected Health+ App

- Bluetooth® 4.2 for automatic data transfer after each measurement
- More on www.microlife.com/connect



MAM: Microlife Average Mode

- 3 measurements in 1 click and average result for superior reliability
- Reliable hypertension diagnosis and therapy control



Clinically Validated

Highest measuring accuracy with clinical A/A grade according to the British and Irish Hypertension Society (BIHS).

12+ years

Validated and recommended for children. Also validated and recommended for patients with diabetes, renal disease, pregnancy / pre-eclampsia and obesity.



The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Microlife Corp. is under license. Other trademarks and trade names are those of their respective owners.

Microlife AG
Epenstrasse 139
9443 Widnau / Switzerland
www.microlife.com

CE0044

www.microlife.com/b6-connect



5 year
guarantee