



Smoothbore catheter mounts

Clear and crush-resistant tubing



Airway Management ▪ Patient Connections

Smoothbore catheter mounts


A full range of Smoothbore catheter mounts for patient connection providing good flow characteristics, lower tubing compliance and clear, crush-resistant tubing. Available with a variety of patient elbows and port cap options.

A guide to our ports and caps

For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.

Luer ports

For secure attachment of gas sampling lines with luer lock fittings




7.6mm port and flip top cap

For insertion of all standard 7.6mm probes or for suctioning



Double flip top cap with seal

Bronchoscopy and suction ports in one patient connection





Code	5120000 (S*)	5180000 (S*)	5250000 (S*)
Box Qty.	50 (85*)	50 (85*)	50 (50*)
Length	120mm	180mm	250mm
Compressible Volume	30ml	37ml	51ml

Sterile option available


(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5120000S (sterile box quantity is shown in brackets).

		
Code	5180104 (S*)	
Box Qty.	50 (85*)	
Compressible Volume	25ml	

	
Code	5180197 (S*)
Box Qty.	50 (85*)
Compressible Volume	35ml

	
Code	5180114 (S*)
Box Qty.	50 (85*)
Compressible Volume	36ml

	
Code	5180193 (S*)
Box Qty.	50 (85*)
Compressible Volume	35ml

	
Code	5180199 (S*)
Box Qty.	50 (85*)
Compressible Volume	37ml



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK

T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.

IS6.14 • Issue 3 07.19