



Patient Connections




Double swivel elbows

A range of double swivel elbows providing added mobility at the patient connection, with a variety of ports for all clinical applications. For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.

Please see page 4 for our guide to ports and caps.

▶ Video available at www.intersurgical.com

	
Code	1897000 (S*)
Box Qty.	35 (15*)
Compressible Volume	12ml

	
Code	1899000 (S*)
Box Qty.	30 (15*)
Compressible Volume	12ml

	
Code	1898000 (S*)
Box Qty.	35 (15*)
Compressible Volume	12ml





Fixed elbows


A range of fixed elbow connections with a variety of ports for a number of clinical applications. For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.

Please see page 4 for our guide to ports and caps.


▶ Video available at www.intersurgical.com


	
Code	1997000 (S*)
Box Qty.	30 (15*)
Compressible Volume	10ml


	
Code	1993000
Box Qty.	30
Compressible Volume	10ml


	
Code	2714000 (S*)
Box Qty.	25 (15*)
Compressible Volume	11ml




	
Code	1883000
Box Qty.	70


	
Code	1992000
Box Qty.	25


	
Code	1991000
Box Qty.	20

	
Code	1893000
Box Qty.	25

	
Code	1891000
Box Qty.	25

	
Code	1892000
Box Qty.	25

	
Code	1995000
Box Qty.	25

	
Code	1990000
Box Qty.	25

◻ Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S (sterile box quantity is shown in brackets).


Flexible catheter mounts


A range of flexible catheter mounts available with a variety of elbow connections and ports for all applications. For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.


Please see page 4 for our guide to ports and caps.


▶ Video available at www.intersurgical.com





	
Code	3511000
Box Qty.	50
Compressible Volume	44ml


Nasal catheter mount 	
Code	3517000
Box Qty.	50
Compressible Volume	43ml


	
Code	3501000 (S*)
Box Qty.	50 (85*)
Compressible Volume	45ml

	
Code	3508000
Box Qty.	50
Compressible Volume	49ml

	
Code	3515000
Box Qty.	40
Compressible Volume	50ml

	
Code	3506000
Box Qty.	50
Compressible Volume	48ml

	
Code	3516000 (S*)
Box Qty.	50 (85*)
Compressible Volume	47ml

	
Covers 7.6mm and 9.5mm applications	
Code	3505000 (S*)
Box Qty.	50 (85*)
Compressible Volume	48ml

Ⓢ Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S (sterile box quantity is shown in brackets).

Superset™ catheter mounts

Superset catheter mounts have a concertina design to allow for complete positioning control.


The design allows the adjustment of compressible volume and the adjustment of length from 70mm to 150mm, excluding connectors.


The range is available with fixed and double swivel elbows and a variety of ports for all applications. For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.


Please see page 4 for the guide to ports and caps.





Video available at www.intersurgical.com


	
Code	3504000 (S*)
Box Qty.	75 (75*)
Compressible Volume	22-43ml


	
Code	3502000 (S*)
Box Qty.	75 (40*)
Compressible Volume	23-44ml


	
Code	3514000 (S*)
Box Qty.	75 (40*)
Compressible Volume	26-47ml

	
Code	3512000
Box Qty.	75
Compressible Volume	27-47ml

	
Code	3520000 (S*)
Box Qty.	75 (40*)
Compressible Volume	23-47ml

	
Code	3521000 (S*)
Box Qty.	75 (40*)
Compressible Volume	26-50ml

Paediatric/neonatal 	
Code	3535000
Box Qty.	15
Compressible Volume	13-16ml

Neonatal 	
Code	3534000
Box Qty.	25
Compressible Volume	10-13ml

 Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S (sterile box quantity is shown in brackets).


Smoothbore catheter mounts


A range of Smoothbore catheter mounts providing good flow characteristics, lower tubing compliance and clear, crush-resistant tubing. Available in a variety of lengths with double swivel elbows and double flip top caps. For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.


Please see page 4 for the guide to ports and caps





Video available at www.intersurgical.com


	
Code	5120000 (S*)
Box Qty.	50 (85*)
Compressible Volume	30ml


	
Code	5180000 (S*)
Box Qty.	50 (85*)
Compressible Volume	37ml


	
Code	5180104 (S*) NEW
Box Qty.	50 (85*)
Compressible Volume	25ml

	
Code	5180114 (S*) NEW
Box Qty.	50 (85*)
Compressible Volume	36ml

	
Code	5180193 (S*) NEW
Box Qty.	50 (85*)
Compressible Volume	35ml

	
Code	5180197 (S*) NEW
Box Qty.	50 (85*)
Compressible Volume	35ml

	
Code	5180199 (S*) NEW
Box Qty.	50 (85*)
Compressible Volume	35ml


	
Code	5250000 (S*)
Box Qty.	50 (50*)
Compressible Volume	51ml


Clear patient connections – reprocessable


A range of reprocessable clear, latex-free connectors with added grip feature to help with the 'push & twist' method to ensure a good gas-tight fit.

The reusable range of polysulfone connectors are designed for use with silicone and Hytrel® tubing. The transparent material allows excellent visibility through the connector. Good dimensional stability over a broad temperature range makes them ideal for reprocessing by steam sterilisation.



	
Code	9114000
Box Qty.	5

	
Code	9115000
Box Qty.	5

	
Code	9120000
Box Qty.	5

 Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S (sterile box quantity is shown in brackets).

Hytrel is a registered trademark of DuPont.