

Single-Use Fibre-Optic Systems

Rüsch Polaris Single-Use Laryngoscope Blade



- Metal Blade construction with the feel of a reusable blade
- Low-Profile blade design to improve visualisation and reduce the risk of dental trauma
- 4.7 mm enclosed Fibre-Optic Cable for greater light transmission
- Colour-coded packaging for easy identification of the correct size

Rüsch Polaris Single-Use Laryngoscope Blade

REF.	DESCRIPTION	OVERALL LENGTH (MM)	DISTAL WIDTH (MM)	COLOUR CODE	QTY
Macintosh blades					
4150100	Mac 0	85.0	10.7	Peach	10
4150110	Mac 1	94.5	9.9	White	
4150120	Mac 2	111.0	13.0	Blue	
4150130	Mac 3	133.0	14.0	Yellow	
4150140	Mac 4	149.00	13.0	Pink	
4150150	Mac 5	159.00	14.5	Dark Yellow	
Miller blades					
4150003	Miller 000	60.0	7.0	Light Green	10
4150002	Miller 00	66.5	12.8	Beige	
4150001	Miller 0	78.0	12.8	Violet	
4150010	Miller 1	103.0	12.8	Orange	
4150020	Miller 2	154.0	13.5	Grey	
4150030	Miller 3	195.0	13.5	Green	
4150040	Miller 4	205.0	19.0	Lavender	

Focused, fibre-optic light transmission to aid visualisation

Performance of single-use with the durability of a reusable blade

Low-Profile Blade Design: Reduced posterior flange height designed to improve visualisation and reduce the risk of dental trauma