

ThinPrep® PreservCyt Vials (Non-GYN)



Item No: ASY-14753

Item Desc: 100-PACK VIAL, NON GYN, INTERNATIONAL

100 vial package, pre-filled with 20ml

ThinPrep® Non-GYN filters



Item No: 70205-001

Item Desc: FILTER, THINPREP NON-GYNECOLOGICAL 100-PACK

Box of 100 filters,
(preferred filter for urine samples)

ThinPrep® UroCyt filters



Item No: 70472-001

Item Desc: FILTER, THINPREP UROCYTE 100-PACK

Box of 100 filters

ThinPrep® Microscope Slides (Non-GYN)



Item No: 70372-001

Item Desc: THINPREP MICROSCOPE SLIDES, NON-GYN, 100-PACK

Box of 100 pc

ThinPrep® Imaging System Slide kit



Item No: 70825-001

Item Desc: SLIDE KIT, THINPREP IMAGING, SMS

Box of 5 x 100 pc, 500-pack

ThinPrep® Microscope, Arcless Slides (Non-GYN)



Item No: 70126-002

Item Desc: 72-SLIDE PACKAGE, THINPREP MICROSCOPE, ARCLSS NON GYN CYTOLOGY

Box of 72

CytoLyt® Solution Bottles



Item No: 70408-002

Item Desc: CYTOLYT 32OZ 4-PACK INTERNATIONAL

4 bottle pack (4 x 946ml)
Transport medium and wash solution

PreservCyt® Solution Bottles



Item No: 70406-002

Item Desc: PRESERVCYT 32OZ 4-PACK INTERNATIONAL

4 bottle pack (4 x 946ml)
Fixative solution

Dispenser Pump



Item No: 50705-001

Item Desc: DISPENSER PUMP 1-PACK

Use for PreservCyt or CytoLyt
Dispensing 30ml

CytoLyt® Solution Collection Cups



Item No: 0236050

Item Desc: CYTOLYT CUPS, 50-PACK

Box of 50 cups. Pre-filled with 30ml
Use for respiratory samples

CytoLyt® Solution Centrifuge Tubes



Item No: 0236080

Item Desc: CYTOLYT TUBES, 80-PACK

Box of 80 centrifuge tubes. Pre-filled 30ml
collection tube for Non-GYN samples

Stackable ThinPrep® Vial Tray



Item No: 70264-001

Item Desc: SAMPLE VIAL TRAY 4-PACK ASSAY

Box of 4 trays
Each tray can hold 40 vials

ThinPrep® Stains



Item No: 70897-002

Item Desc: THINPREP STAIN SET, INTERNATIONAL

Kit containing 5 stains:

- 4 liters Rinse Solution (70779-001)
- 4 liters Nuclear Stain (70780-001)
- 4 liters Orange G Solution (70781-002)
- 4 liters EA Solution (70782-002)
- 4 liters Bluing Solution (70793-001)

ThinPrep® Stains



Item No: 70779-001

**Item Desc: THINPREP RINSE SOLUTION,
4L PACKAGED PRODUCT, INTERNATIONAL**

ThinPrep® Stains



Item No: 70780-001

**Item Desc: THINPREP NUCLEAR STAIN,
4L PACKAGED PRODUCT, INTERNATIONAL**

ThinPrep® Stains



Item No: 70781-002

**Item Desc: THINPREP ORANGE G SOLUTION
4L PACKAGED PRODUCT, INTERNATIONAL**

ThinPrep® Stains



Item No: 70782-002

**Item Desc: THINPREP EA SOLUTION 4L
PACKAGED PRODUCT, INTERNATIONAL**

ThinPrep® Stains



Item No: 70793-001

**Item Desc: THINPREP BLUING SOLUTION
4L PACKAGED PRODUCT, INTERNATIONAL**