



GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: Surgical Probe 16cm, 1.5mm (KNAPPSOND)



Glorious Art nr.	54.0162.16
Finish:	Satin
Net weight:	0.002 gram (approximate)
Material:	Stainless steel AISI 410
Sterility:	Non-Sterile (Re-Usable)
Exposure Temperature:	134°C (274°F)
Exposure Time:	3 Minutes
CE Classification:	Class-I
Storage:	16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)
life-span:	20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.
Disposal:	No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material



P.O. Box 63, Naverland 11, DK 2600 Glostrup, Denmark
Tel: +45 43 62 26 93, Fax: +45 43 62 26 11 , Email: glorious@glorious-enterprises.dk