

# BD SmartSite™ Extension Set

**SKU/REF** MFX2239E  
**SKU/REF Name** BD SmartSite™ 3-Way Extension Set - 3 BD SmartSite™ - 3 back-check valves - 12 cm



## Specifications

Tube Length (cm)/ Total Length (cm)	7/12
Priming Volume (mL)	0.6
Tube Inner Diameter (mm)	1
Ways	3-Way
Clamp/ Back-Check Valve	3/3
Luer Lock	Rotating Male
Needle Free Connector	BD SmartSite™
Gravity Flow Rate	85 mL/min
Pressure Rating	Max. 2.5 bar
PVC	Yes
PU	Yes
DEHP	DEHP-Free
Shelf Life	5 Years
Sterilization Method	Ethylene Oxide
Country Of Manufacture	Malta
Classification	Ila
Quantity per Box/Case	50/150

CE 0481

Becton Dickinson Sverige AB, Marieviksgatan 25, 11743 Stockholm

[bdsweden@bd.com](mailto:bdsweden@bd.com)

[bd.com](https://www.bd.com)

BD, BD-logo och SmartSite är varumärken som tillhör Becton, Dickinson and Company eller dess dotterbolag  
© 2022 BD. Alla rättigheter förbehållna. BD-59783 March 2022

