



### **Product description**

Order number	62.553.542
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, yellow, cap assembled, with print, label/print: white/blue, with graduations, sterile, 50 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,  
with print**

**Page 2**

**Size**

Diameter 17 mm  
Length of product 120 mm

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow

**Purity & certification**

Satisfies the requirement sealed at 0.6 bar/10 min  
Product category In vitro diagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Pyrogen-free yes  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 24000  
Depth of case 391 mm  
Width of case 219 mm  
Height of case 429 mm  
Case volume 0.0367 cbm  
Weight of product 0.0085 kg  
Weight of case 4.26 kg  
EAN of inner box 4038917113007  
EAN case 4038917013871

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-10-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**