

microlife®

A partner for people. For life.



Key features



3 measurements in 1 click and average result for superior reliability



Irregular Heartbeat (IHB) detection for early alerts to possible heart irregularities



For a accurate and extra comfortable measurement



Easy to read classification of blood pressure result



Micro-USB cable included to transfer your data from the BP B3 to your PC and the included software BP Analyzer+



99 data memory and average-display of all stored readings



Cuff fit check

BP_{B3} Comfort PC

Blood Pressure Monitor

- Technology Comfort+ for measurement during inflation
- MAM: 3 measurements in 1 click & average result
- Comfortable cuff for most arm sizes (22-42 cm)
- Cuff fit check
- Date and time



2 users

- with 99 data memory each
- average display for all stored readings



Soft Cuff: 22-42 cm

- Conical shape for better fit and optimal accuracy
- Washable cuff cover for longlasting, hygienic use



Clinically Validated BIHS Protocol A/A

Highest measuring accuracy with clinical A/A grade according to the British Hypertension Society (BIHS).

12+
years



Validated and recommended for children. Also validated and recommended for patients with diabetes, renal disease, pregnancy/pre-eclampsia and obesity.



5 year
guarantee