

Isolation gown

distributed by Mölnlycke®



- Impervious material, nonwoven laminated with polyethylene
- Three-threads overlock for sleeve and body-sleeve joints
- Non-sterile

Product classification CE mark Directive 93/42/EEC

Material 18gsm polypropylene nonwoven laminated with 22gsm colored polyethylene

Material hydrohead cmH₂O >100

Colour Blue

Country of origin China

Order information

Ref. No.	Sterile codes	Size	Length (cm)	Width (cm)	Disp/Trp
31102C	Isolation gown	L	115	137	10/100
31103C	Isolation gown	XL	120	137	10/100