

Methionine¹⁰⁰

ACBS ☐ GMS ☐ PBS ☐

Description

Food for special medical purposes.

A powdered methionine amino acid supplement on a carbohydrate base.

Indications

For the dietary management of inborn errors of amino acid metabolism from birth.

Dosage and Administration

To be determined by the clinician or dietitian and is dependent on the age, body weight and medical condition of the patient.

Preparation Guidelines

To be determined by the clinician or dietitian.

May be used in low protein diets in a variety of ways:

- Mixed in with protein substitute
- Added to modular feeds.
- Incorporated into other permitted food or drink.

Methionine100 (4g) = 100mg L- Methionine

IMPORTANT NOTICE

Must only be given to patients with proven inborn errors of amino acid metabolism under strict medical supervision.

Suitable from birth.

Children under 12 months: when made into solution or added to infant formula, use immediately and discard any remaining within 1 hour.

Not for use as a sole source of nutrition.

For enteral use only.

Nutritional Information

		Per 100g	Per 4g sachet		Per 100g	Per 4g sachet
Energy	kJ	1610	64	Protein Equivalent	g	2.2
		379	15	L-Methionine	mg	2500
Fat	g	0	0	Salt	g	<0.05
		0	0			<0.002
of which Saturates	g	0	0			
Carbohydrate	g	92.5	3.7			
of which Sugars	g	8.8	0.4			

Storage

Store in a cool, dry place.

Sachets are designed for single use and should be used once opened.

Pack Size/Weight

30 x 4g sachets = 120g

Flavour

Unflavoured

