



GE Healthcare

BiliSoft™ Disposable Pad Covers and Pad Nests

Intensive therapy, as easy as wrapping a baby in a blanket

The BiliSoft Pad Covers and Pad Nests are designed for use with the BiliSoft Phototherapy Systems, which are used to treat both premature and full-term infants for hyperbilirubinemia.

Features

- Ultra-soft, cushioned, disposable BiliSoft covers are made of flame-retardant, skin-friendly fabric
- Available in small and large sizes to match the size of the light pad
- Materials used are engineered to allow therapy light to pass through the fabric
- Covers are x-ray compatible

**Our equipment.
Our accessories.
Simply compatible.**

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step



gehealthcare.co.uk/products/accessories-and-supplies

PRODUCT SPECIFICATIONS

Product use

BiliSoft Pad covers are gentle on the baby's skin and provide a soft cushioned covering for all BiliSoft products.

Materials

Top Layer:	Non-woven polypropylene, no dye
Inner Padding:	Polyester fiber batting, no dye
Back Panel:	Polyester Mesh, no dye
Stitching:	Polyester thread

Materials composition

Not made with any BPA, DEHP, natural rubber latex, PVC or other phthalates.

Operating environment

Temperature:	5° to 40°C (41 to 104°F)
Humidity:	up to 95% RH (non-condensing)

Storage

Temperature:	-25° to 70°C (-13° to 158°F)
Humidity:	up to 90% RH (non-condensing)

Disposal considerations

The BiliSoft disposable covers are for single-patient-use only. Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.co.uk.

Packaging materials

Primary: Polybag – Polyethylene (PE)

Regulatory standards

Pad Covers Meet: ISO 10993-5 and ISO 10993-10 for Biocompatibility Requirements and 16 CFR Part 1632.6 for Flammability Requirements.

Compliance

CE Mark, regulation (EU) 2017/745, REACH regulation 1907/2006 compliant, including amendments.

PART LIST

Part Number	Description
2105422-001	BiliSoft 2.0 Disposable Pad Covers, Flat, Small, 20pk
2105418-001	BiliSoft 2.0 Disposable Pad Covers, Flat, Small, 50pk
2107201-001	BiliSoft Disposable Pad Covers, Flat, Small, 20pk
2107200-001	BiliSoft Disposable Pad Covers, Flat, Small, 50pk
2105423-001	BiliSoft 2.0* Disposable Pad Covers, Flat, Large, 20pk
2105419-001	BiliSoft 2.0* Disposable Pad Covers, Flat, Large, 50pk
2105420-001	BiliSoft 2.0 Disposable Nest Pad Covers, Small, 15pk
2107202-001	BiliSoft Disposable Nest Pad Covers, Small, 15pk
2105421-001	BiliSoft 2.0* Disposable Nest Pad Covers, Large, 15pk

*Note: the BiliSoft 2.0 Disposable Pad (flat and nest) in size large also fits the BiliSoft Phototherapy System (legacy) large pad.

Distributed by:
GE Healthcare



Manufactured by:
Lumitex Medical Devices Inc.
8443 Dow Circle
Strongsville, OH 44136