

microlife®

A partner for people. For life.

Small, light and
PORTABLE

Simple and completely
pain-free measurement



OXY 500 BT

Fingertip Pulse Oximeter

- Determining oxygen saturation (SpO₂) and heart rate (pulse)
- Comfortable measurement
- Bluetooth® technology for easy monitoring

Key features



Medical grade accuracy

98 75



OLED display with SpO₂ and pulse display and bar graph display



Simple and pain-free measurement



Easy to read - Display with 6 available views



Instant monitoring on your smartphone with the Microlife Connected Health+ app



Instant tracking on your smartphone



Ideal for use at home and on the move



Automatic switch-off



Adjustable display brightness



Risk reminder



CE 0482

2 year
guarantee