


BLOOD COLLECTION



Sol-Care™ Safety Blood Collection Needle



Safety-engineered device intended for venous blood collection. The safety mechanism can be activated immediately after blood draw, helping to protect you against needlestick injury.

It is also offered with a pre-attached holder.

Sol-Care™

● Safety Blood Collection Needle



CE 0123

Sol-Care™ Safety Blood Collection Needle with Luer Adapter

Sol-Care™ Safety Blood Collection Needle with Pre-attached Holder



Sol-Care™ Safety Blood Collection Needle

with Luer Adapter

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Tubing Length (nches)	Tubing Length (mm)	Box Qty	Case Qty
110201030035	25 G	0.5 mm	3/4"	19 mm	12"	305 mm	50	200
110201030033	23 G	0.6 mm	3/4"	19 mm	12"	305 mm	50	200
110201030031	21 G	0.8 mm	3/4"	19 mm	12"	305 mm	50	200
110201030034	25 G	0.5 mm	3/4"	19 mm	7"	178 mm	50	200
110201030032	23 G	0.6 mm	3/4"	19 mm	7"	178 mm	50	200
110201030030	21 G	0.8 mm	3/4"	19 mm	7"	178 mm	50	200



Sol-Care™ Safety Blood Collection Needle

with pre-attached holder

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Tubing Length (nches)	Tubing Length (mm)	Box Qty	Case Qty
110201030010	25 G	0.5 mm	3/4"	19 mm	12"	305 mm	25	200
110201030006	23 G	0.6 mm	3/4"	19 mm	12"	305 mm	25	200
110201030002	21 G	0.8 mm	3/4"	19 mm	12"	305 mm	25	200
110201030012	25 G	0.5 mm	3/4"	19 mm	7"	178 mm	25	200
110201030008	23 G	0.6 mm	3/4"	19 mm	7"	178 mm	25	200
110201030004	21 G	0.8 mm	3/4"	19 mm	7"	178 mm	25	200