


BLOOD COLLECTION



Sol-M™ Multi-sample Needle **Sol-Care™ Safety Multi-sample Needle**



Low-angle bevel design
Allows easy penetration.

Color-coded caps and hubs
Allow easy identification of needle gauge size.

Not made with natural rubber latex.

Sol - M™

Multi-sample Needle

Sol-M™ Multi-Sample Needle makes it possible to collect several samples with a single vein puncture.

- Marker on needle cap allows for bevel-up identification



Sol-M™ Multi-sample Needle

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Box Qty	Case Qty
MSN2310	23 G	0.6 mm	1"	25 mm	100	1000
MSN2210	22 G	0.7 mm	1"	25 mm	100	1000
MSN2215	22 G	0.7 mm	1 1/2"	38 mm	100	1000
MSN2110	21 G	0.8 mm	1"	25 mm	100	1000
MSN2115	21 G	0.8 mm	1 1/2"	38 mm	100	1000
MSN2010	20 G	0.9 mm	1"	25 mm	100	1000
MSN2015	20 G	0.9 mm	1 1/2"	38 mm	100	1000
VMSN1810	18 G	1.2 mm	1"	25 mm	100	1000
VMSN1815	18 G	1.2 mm	1 1/2"	38 mm	100	1000

18G needle has been designed for veterinary use; contains a component made of natural rubber latex.

Sol-Care™

Safety Multi-sample Needle

Safety-engineered tool for multi-sample venous blood collection. The device is designed with a safety mechanism covering the needle after use. In the activated position, the needle cover guards against accidental needlestick during normal handling and disposal.

It is also offered with a pre-attached holder.

Sol-Care™ Safety Multi-sample Needle

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Box Qty	Case Qty
SMSN22125	22 G	0.7 mm	1 1/4"	32 mm	48	480
SMSN2215	22 G	0.7 mm	1 1/2"	38 mm	48	480
SMSN21125	21 G	0.8 mm	1 1/4"	32 mm	48	480
SMSN2115	21 G	0.8 mm	1 1/2"	38 mm	48	480
SMSN20125	20 G	0.9 mm	1 1/4"	32 mm	48	480
SMSN2015	20 G	0.9 mm	1 1/2"	38 mm	48	480

Sol-Care™ Safety Multi-sample Needle with Pre-attached Holder

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (nches)	Needle Length (mm)	Box Qty	Case Qty
SMSNH22125	22 G	0.7 mm	1 1/4"	32 mm	48	480
SMSNH2215	22 G	0.7 mm	1 1/2"	38 mm	48	480
SMSNH21125	21 G	0.8 mm	1 1/4"	32 mm	48	480
SMSNH2115	21 G	0.8 mm	1 1/2"	38 mm	48	480
SMSNH20125	20 G	0.9 mm	1 1/4"	32 mm	48	480
SMSNH2015	20 G	0.9 mm	1 1/2"	38 mm	48	480



MULTI-SAMPLE NEEDLE

