



Produced in the EU for:



Important Notice

Must only be given to patients with proven PKU (phenylketonuria) under strict medical supervision. Suitable from 3 years of age. Not for use as a sole source of nutrition. Diet must be supplemented with natural protein, water and other nutrients in prescribed quantities to supply phenylalanine, fluid and general nutritional requirements. For enteral use only.


Storage

Store in a cool, dry place.

Pack Size

6 x 77 tablets.

A Food for Special Medical Purposes



PKU Easy Tablets

PKU Easy Tablets are a phenylalanine free protein substitute tablet containing essential and non essential amino acids. PKU Easy Tablets are suitable to be used in the dietary management of PKU (phenylketonuria) patients aged 3 years and above.

Dosage PKU Easy Tablets

Dosage to be determined by the clinician or dietitian and is dependent on the age, body weight and medical condition of the patient. PKU Easy Tablets are coated (slow release) tablets which are to be swallowed with liquid. 11 tablets = 10 gram of protein.

Nutritional information

| Nutritional Information | Unit | 100g |
|-------------------------|------|------|
| Energy | kJ | 1659 |
| | | Kcal |
| Protein Equivalent | g | 73 |
| Amino Acids | g | 88 |
| Carbohydrate | g | 12.5 |
| of which sugars | g | 0 |
| Fat | g | 0 |
| Of which saturated | g | 0 |
| Fiber | g | 0 |

| Amino Acids | Unit | 100g |
|-----------------|------|------|
| L-Alanine | g | 3.2 |
| L-Arginine | g | 5.1 |
| L-Aspartic Acid | g | 8.0 |
| L-Carnitine | g | 0.08 |
| L-Cystine | g | 2.1 |
| L-Glutamine | g | 6.2 |
| Glycine | g | 8.0 |
| L-Histidine | g | 3.2 |
| L-Isoleucine | g | 5.5 |
| L-Leucine | g | 8.6 |
| L-Lysine | g | 5.7 |
| L-Methionine | g | 1.5 |
| L-Phenylalanine | g | - |
| L-Proline | g | 5.7 |
| L-Serine | g | 3.5 |
| Taurine | g | 0.12 |
| L-Threonine | g | 5.5 |
| L-Tryptophan | g | 1.7 |
| L-Tyrosine | g | 8.1 |
| L-Valine | g | 6.3 |