



Because **OUTCOMES** matter  
**BEYOND** the **NICU.**



# fabian™ Ventilators

product catalog

Vyaire Medical Inc. is a global company focused exclusively on supporting breathing through every stage of life. Established from legacy brands with a 65-year history of pioneering breathing technology, the company's portfolio of integrated solutions is designed to enable, enhance and extend lives.

At Vyaire, we pledge to advance innovation and evolve what's possible to ensure every breath is taken to its fullest. Headquartered in suburban Chicago, Vyaire is recognized, trusted and preferred by specialists throughout the respiratory community worldwide.

# Table of Contents

List of Abbreviations .....	4
-----------------------------	---

## Devices

fabian™ HFOi .....	5
fabian™ HFOi Light .....	5
fabian™ +nCPAP evolution .....	6
fabian™ Therapy evolution .....	6

## Software Options

PRICO-NIV Closed Loop Control FiO <sub>2</sub> _SpO <sub>2</sub> .....	7
PRICO-CNV Closed Loop Control FiO <sub>2</sub> _SpO <sub>2</sub> .....	7
PRICO-CNV Closed Loop Control FiO <sub>2</sub> .....	7
PRICO-HFO Closed Loop Control FiO <sub>2</sub> _SpO <sub>2</sub> .....	7
NIV Trigger for fabian™ Therapy evolution .....	8
NIV Trigger for fabian™ +nCPAP & HFO .....	8
Lung Recruitment Software for HFO-Mode.....	8
FOT Software Option License Key.....	8
PDMS Licence Key .....	8

## Accessories

Upgrade Kits.....	9
fabian™ Flow Sensors .....	9
Exhalation Valve Parts.....	10
Testlungs .....	10
Humidification .....	10-12
Dual Limb Ventilation.....	12-13
Single Limb Ventilation.....	13-14
Patient Circuit Spares.....	14
Infant Flow™ nCPAP Generator .....	14
Infant Flow™ LP Headgear.....	15
Infant Flow™ LP Bonnets.....	15
Infant Flow™ LP Nasal Prongs .....	15
Infant Flow™ LP Nasal Mask .....	16
High Flow Nasal Cannula.....	16
SpO <sub>2</sub> Option Masimo .....	17
etCO <sub>2</sub> Option Microstream .....	17-18
etCO <sub>2</sub> Option Lo-Flow.....	18
Aerogen Solo USB Version.....	19
Aerogen PRO Device .....	20-21
Carts.....	21
Air Compressors.....	22-23
Carts Accessories.....	24
Carts Spares.....	24
Extended Warranty Program .....	25
Air and O <sub>2</sub> Sets .....	25-26
Powercords Angled .....	26-27

# List of Abbreviations

<b>AI</b>	fabian™ HFOi
<b>AH</b>	fabian™ HFO Classic
<b>AL</b>	fabian™ HFOi Light
<b>AK</b>	fabian™ HFO Light Classic
<b>AN</b>	fabian™ +nCPAP evolution
<b>AT</b>	fabian™ Therapy evolution
<b>+nCPAP</b>	fabian™ +nCPAP (old version)
<b>Therapy</b>	fabian™ Therapy (old version)



**Part no. Description**

113001

**fabian™ HFOi  
Neonatal & Paediatric HFO Ventilator**

**Battery backup:** Nickel-Metal-Hydride-Battery  
**Backup time:** Conventional mode 2.5 hours, HFO mode 1 hour  
**User interface:** Colour touchscreen  
**Mode:** Time cycled, pressure controlled  
**Ventilation modes:** HFO, IPPV/IMV, SIPPV, SIMV, SIMV+PSV, PSV, CPAP, NCPAP, DUOPAP, O<sub>2</sub> Therapy, Manual breath  
**Included enhancements:** Volume Limit (VL) & Volume Guarantee (VG), data trend 5 days  
**Power supply:** 100-240 V, 50/60 Hz  
**Gas connection type (air and oxygen):** NIST F

Including:

- 1 pc x SW module NCPAP and DUOPAP REF 150212
- 1 pc x SW module O<sub>2</sub> high flow therapy REF 150211
- 1 pc x SW module data trend 5 days REF 150226
- 1 pc x SW module VG & VL REF 150221
- 1 x Accessory kit "invasive" REF 100001.02 containing itself:
  - > 1 pc x Flow sensor cable, reusable REF 151100.01
  - > 1 pc x Flow sensor MPU neonatal REF 151120
  - > 1 pc x Testlung REF 151900
  - > 2 pcs x Infant Flow™ connecting tube REF 153011.02

Please choose a power cord when you place the order!

- 1 pc x Power cord (REF 497180.XX)



111001.01

**fabian™ HFOi Light (without HFO module)  
Neonatal & Paediatric Invasive Unit**

**Battery backup:** Nickel-Metal-Hydride-Battery  
**Backup time:** Conventional mode 2.5 hours  
**User interface:** Colour touchscreen  
**Mode:** Time cycled, pressure controlled  
**Ventilation modes:** IPPV/IMV, SIPPV, SIMV, SIMV+PSV, PSV, CPAP, NCPAP, DUOPAP, O<sub>2</sub> Therapy, Manual breath  
**Included enhancements:** Volume Limit (VL) & Volume Guarantee (VG), data trend 5 days  
**Power supply:** 100-240 V, 50/60 Hz  
**Gas connection type (air and oxygen):** NIST F

Including:

- 1 pc x SW module NCPAP and DUOPAP REF 150212
- 1 pc x SW module O<sub>2</sub> high flow therapy REF 150211
- 1 pc x SW module data trend 5 days REF 150226
- 1 pc x SW module VG & VL REF 150221
- 1 x Accessory kit "invasive" REF 100001.02 containing itself:
  - > 1 pc x Flow sensor cable, reusable REF 151100.01
  - > 1 pc x Flow sensor MPU neonatal REF 151120
  - > 1 pc x Testlung REF 151900
  - > 2 pcs x Infant Flow™ connecting tube 153011.02

Please choose a power cord when you place the order!

- 1 pc x Power cord (REF 497180.XX)



Part no.	Description
----------	-------------

122001

**fabian™ +nCPAP evolution  
Neonatal & Paediatric Invasive Unit**

**Battery backup:** Nickel-Metal-Hydride-Battery  
**Backup time:** 3 hours  
**User interface:** Colour touchscreen  
**Mode:** Time cycled, pressure controlled  
**Ventilation modes:** IPPV/IMV, SIPPV, SIMV, SIMV+PSV, PSV, CPAP, NCPAP, DUOPAP, O<sub>2</sub> Therapy, Manual breath  
**Included enhancements:** Volume Limit (VL) & Volume Guarantee (VG), data trend 5 days  
**Power supply:** 100-240 V, 50/60 Hz  
**Gas connection type (air and oxygen):** NIST F

Including:

- 1 pc x SW module NCPAP and DUOPAP REF 150212
- 1 pc x SW module O<sub>2</sub> high flow therapy REF 150211
- 1 pc x SW module data trend 5 days REF 150226
- 1 pc x SW module VG & VL REF 150221
- 1 x Accessory kit "invasive" REF 100001.02 containing itself:
  - > 1 pc x Flow sensor cable, reusable REF 151100.01
  - > 1 pc x Flow sensor MPU neonatal REF 151120
  - > 1 pc x Testlung REF 151900
  - > 2 pcs x Infant Flow™ connecting tube 153011.02

Please choose a power cord when you place the order!

- 1 pc x Power cord (REF 497180.XX)



121001

**fabian™ Therapy evolution  
NIV Ventilator**

**Battery backup:** Nickel-Metal-Hydride-Battery  
**Backup time:** 3 hours  
**User interface:** Colour touchscreen  
**Mode:** Time cycled, pressure controlled  
**Ventilation modes:** NCPAP, DUOPAP, O<sub>2</sub> Therapy, Manual breath  
**Included enhancements:** Data trend 5 days  
**Power supply:** 100 - 240 V, 50/60Hz  
**Gas connection type (air and oxygen):** NIST F







Including:

- 1 pc x SW module NCPAP and DUOPAP REF 150212
- 1 pc x SW module O<sub>2</sub> high flow therapy REF 150211
- 1 pc x SW module data trend 5 days REF 150226
- 2 pcs x Infant Flow™ connecting tube REF 153001

Please choose a power cord when you place the order!

- 1 pc x Power cord (REF 497180.XX)



Part no.	Description	
150216	<b>PRICO-NIV Closed Loop Control <math>\text{FiO}_2\text{-SpO}_2</math></b> <b>For: AT</b> 4th Edition – including cable LEMO 152406 and SW key 150238	
150217	<b>PRICO-CNV Closed Loop Control <math>\text{FiO}_2\text{-SpO}_2</math></b> <b>For: AK, AL, AN</b> 4th Edition – including cable LEMO 152406 and SW key 150239	
150218	<b>PRICO-HFO Closed Loop Control <math>\text{FiO}_2\text{-SpO}_2</math></b> <b>For: AI, AH</b> 4th Edition – including cable LEMO 152406 and SW key 150240	
150210	<b>PRICO-NIV Closed Loop Control <math>\text{FiO}_2</math></b> <b>For: AT</b> 3rd Edition - incl. Masimo Module USB (REF 152401) and SW key 150238	
150220	<b>PRICO-CNV Closed Loop Control <math>\text{FiO}_2</math></b> <b>For: AK,AL,AN</b> 3rd Edition - incl. Masimo Module USB (REF 152401) and SW key 150239	
150230	<b>PRICO-HFO Closed Loop Control <math>\text{FiO}_2</math></b> <b>For: AI, AH</b> 3rd Edition - incl. Masimo Module USB (REF 152401) and SW key 150240	




Part no.	Description
150213	<p><b>NIV Trigger for fabian™ Therapy evolution</b> <b>For: AT</b></p> <p>This option enables the NIV flow trigger functionality for devices manufactured AFTER April 2017 (SN higher than AT-001920). Must be used with software version &gt;=5.0 Option for trigger in NCPAP/DUOPAP modes with single limb ventilation and apnea detection in nCPAP mode Flow triggering (for Infant Flow™ &amp; Infant Flow™ LP generators) in DUOPAP Pressure triggering (for Medijet generator) in DUOPAP Allowing S-DUOPAP functionalities</p> <p>Including:</p> <ul style="list-style-type: none"><li>- 1 pc x Software licence key REF 150215</li><li>- 1 pc x Flow sensor cable, reusable REF 151100.01</li><li>- 1 box (à 10 pcs) x Flow sensor, SPU, neonatal REF 151111</li></ul>
150223	<p><b>NIV Trigger for fabian™ +nCPAP &amp; HFO</b> <b>For: AI,AH,AK,AL,AN</b></p> <p>Option for NCPAP/DUOPAP modes in single limb ventilation Apnea detection in nCPAP mode Flow triggering (for Infant Flow™ &amp; Infant Flow™ LP generators) in DUOPAP mode Pressure triggering (for Medijet generator) in DUOPAP mode Only for software version &gt;=5.0</p>
150234	<p><b>Lung Recruitment Software for HFO-Mode</b> <b>For: AI,AH</b></p> <p>Function to open the lung by increasing the pressure level cyclically</p>
150236	<p><b>FOT Software Option License Key</b> <b>For: AI,AH</b></p> <p>Option available for GUI Softwares &gt;5.0 Best PEEP determination tool by Forced Oscillation Technique (FOT)</p>
150225	<p><b>PDMS Licence Key</b> <b>For: AI,AH,AK,AL,AN,AT</b></p> <p>For Vuelink, AcuLink and Intellibridge</p>





## Upgrade Kits

Part no.	Description
150214	<p><b>NIV Trigger Upgrade Kit</b> <b>For: AT</b></p> <p>This kit enables the NIV flow trigger functionality for devices manufactured BEFORE April 2017 (SN lower than AT-001920) Must be used with software version &gt;=5.0 Option for trigger in NCPAP/DUOPAP modes with single limb ventilation and apnea detection in nCPAP mode Flow triggering (for Infant Flow™ &amp; Infant Flow™ LP generators) in DUOPAP Pressure triggering (for Medijet generator) in DUOPAP Allowing S-DUOPAP functionalities Including:</p> <ul style="list-style-type: none"> <li>- 1 pc x Software licence key REF 150215</li> <li>- 1 pc x Connector PCB for AT device REF 122222</li> <li>- 1 pc x Main board rev. 2.2.1 REF 121200.22</li> <li>- 1 pc x Flow sensor cable, reusable REF 151100.01</li> <li>- 1 box (à 10 pcs) x Flow sensor, SPU, neonatal REF 151111</li> </ul>



## fabian™ Flow Sensors

Part no.	Description	
151100.01	<p><b>Flow Sensor Cable, Reusable</b> <b>For: AI,AH,AK,AL,AN,AT,+nCPAP</b></p> <p>Cable lenght: 2 m Weight: 82.1g Compatible with all in this pricelist listed flow sensors</p>	
151111	<p><b>Flow Sensor, SPU, Neonatal 10/box</b> <b>For: AI,AH,AK,AL,AN,AT,+nCPAP</b></p> <p>Weight: 10.5g Compatible with flow sensor cable REF 151100.01</p>	
151120	<p><b>Flow Sensor, MPU, Neonatal</b> <b>For: AI,AH,AK,AL,AN,AT,+nCPAP</b></p> <p>Autoclavable at 134° C Individually packed Weight: 10.9g Compatible with flow sensor cable REF 151100.01</p>	
151121	<p><b>Flow Sensor, MPU, Paediatric</b> <b>For: AI,AH,AK,AL,AN,AT</b></p> <p>Autoclavable at 134° C Individually packed Compatible with flow sensor cable REF 151100.01</p>	


## Exhalation Valve Parts

Part no.	Description	
151501.01	<b>Exhalation Valve Membrane, MPU</b> <b>For: AI,AH,AK,AL,AN,+nCPAP</b> Autoclavable at 134° C	
151520.01	<b>Membrane Holder, 22OD/15ID, MPU, PC</b> <b>For: AI,AH,AK,AL,AN,+nCPAP</b> Made of polycarbonate (PC) Autoclavable at 134°C for 30 cycles maximum.	






## Testlungs

Part no.	Description	
151900	<b>Testlung, 1 Bellows, Neonatal</b> <b>For: AI,AH,AK,AL,AN,+nCPAP</b> Material: Silicone Autoclavable at 134° C	
151901	<b>Testlung, 2 Bellows, Paediatric</b> <b>For: AI,AH,AK,AL,AN,+nCPAP</b> Test lung for neonatal ventilators, with 2 bellows, to simulate one lung ventilation and ET tube leakage	




## Humidification

Part no.	Description	
152120.01EU	<b>Humidifier, WILamed, AIRcon G2, 220V/50H</b> <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b> Colour LCD respiratory gas humidifier with dual servo control for airway temperature Including: - 1 pc x Temperature probe - 1 pc x Heater wire adapter (i+e, Wilamed type) - 1 pc x Fixation for bracket - 1 pc x User manual - 1 pc x Short instructions - 1 pc x Power cord, type C (for Europe) Please choose the correct power cord when type C is not suitable! - 1 pc x Power Cords Angled (REF 497090.XX)	



## Humidification (cont.)

Part no.	Description	Image
152120.01UK	<p><b>Humidifier WILamed, AIRcon G2, 220V/50Hz</b>  <b>For: AI,AH,AK,ALAN,AT,+nCPAP,Therapy</b></p> <p>Colour LCD respiratory gas humidifier with dual servo control for airway temperature</p> <p>Including:</p> <ul style="list-style-type: none"> <li>- 1 pc x Temperature probe</li> <li>- 1 pc x Heater wire adapter (i+e, Wilamed type)</li> <li>- 1 pc x Fixation for bracket</li> <li>- 1 pc x User manual</li> <li>- 1 pc x Short instructions</li> <li>- 1 pc x Power cord, type G (for UK)</li> </ul> <p>If another power cord is required, please order REF 152120.01EU</p>	
152121	<p><b>Temperature Probe for Aircon 150 cm</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Dual way</p>	
152124	<p><b>Dual Heated Wire Adapter to F&amp;P</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For WILAMED humidifier, model AIRcon, compatible with ALTERA and AirLife™ breathing circuits</p> <p>Clover leaf type connectors</p>	
152141	<p><b>Adapter Cable Dual Heater Wires LEMO</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Inspiratory + expiratory heater wire adapter cable</p> <p>To be used with humidifier VH3000</p> <p>Can not be used with WILamed AIRcon humidifier</p> <p>For Multiple Patient Use circuits</p> <p>LEMO connector</p>	
152142	<p><b>Temperature Probe 150cm for VH3000 MPU</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>To be used with humidifier VH3000 (REF 152140)</p> <p>Length: 150 cm</p>	



## Humidification (cont.)

Part no.	Description	
152143	<p><b>Adapter Cable Dual Heater Wires F&amp;P</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Inspiratory + expiratory heater wire adapter cable                      To be used with humidifier VH3000                      F&amp;P connector</p>	
152190	<p><b>Neo Humidifier Chamber MPU</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>To be used with humidifier VH3000                      Autoclavable</p>	
AH290	<p><b>Humidification chamber, automatic feed, SPU 10/Box</b></p> <p>Single use</p>	





## Dual Limb Ventilation

Part no.	Description	
154001	<p><b>Neo Dual Heated Circuit MPU LEMO</b>  <b>For: AI,AH,AK,AL,AN,+nCPAP</b></p> <p>Silicone tubing circuit with LEMO dual heated plugs                      To be used with humidifier VH3000 REF 152140 and LEMO adapters REF 152141                      Without humidification chamber                      Autoclavable</p>	
154308.02	<p><b>Neo Tubing Set, Dual Heated, SPU 10/Box</b>  <b>For: AI,AH,AK,AL,AN,+nCPAP</b></p> <p>Made by ALTERA                      Including:</p> <ul style="list-style-type: none"> <li>- 1 pc x Proximal pressure line</li> <li>- 1 pc x Incubator extension tube</li> <li>- 1 pc x Infant Flow™ connecting tube</li> <li>- 1 pc x Various adapter to convert to NIV ventilation</li> <li>- 1 pc x Autofeed AL-0005 humidifier chamber (NEW)</li> </ul>	


## Dual Limb Ventilation (cont.)

Part no.	Description	
AH165	<p><b>Infant dual-limb, dual-heated circuit, SPU 10/box</b>  <b>For: AI,AH,AK,AL,AN,+nCPAP</b>                      4' high flow circuit (&gt; 4 L/min), pressure line, dry gas line and adapters</p>	
AH265	<p><b>Infant dual-limb, dual-heated circuit with chamber, SPU 10/box</b>  <b>For: AI,AH,AK,AL,AN,+nCPAP</b>                      4' high flow circuit with chamber (&gt; 4 L/min), pressure line, dry gas line and adapters</p>	
154902	<p><b>iNO Kit for fabian™ HFO SPU 10/Box</b>  <b>For: AI,AH,AK,AL,AN,AT</b>                      With Luer-lock adapters, one-way valve and 22F-22F straight connector                      For connection to e.g. a NOxBOX system</p>	




## Single Limb Ventilation

Part no.	Description	
153011.02	<p><b>Infant Flow™ Tube, 60cm, SPU, 10/Box</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p>	
AH032	<p><b>AirLife™ Infant single-limb circuit, SPU 10/box</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b>                      Heated 4' continuous flow circuit (&gt; 4 L/min), pressure line, dry gas line and adapters</p>	
AH132	<p><b>AirLife™ Infant single-limb circuit with chamber, SPU 10/box</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b>                      Heated 4' continuous flow circuit with chamber (&gt; 4 L/min), pressure line, dry gas line and adapters</p>	
153800	<p><b>Transport nCPAP Tubing Set SPU, 10/Box</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b>                      Box consisting of:                      - 10 pcs x Transport single use Infant Flow™/INSPIRE connecting tube                      - 10 pcs x Hose between the unit and flow generator</p>	
156001	<p><b>nCPAP Adapter Kit SPU, 15/box</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b>                      For Infant Flow™ and Infant Flow™ LP generators</p>	


## Single Limb Ventilation (cont.)

Part no.	Description	
156002	nCPAP Adapter Kit SPU, 15/box For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy	

## Patient Circuit Spares

Part no.	Description	
155200	Connector 15OD/9-11OD SPU, 10/Box For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy	
156206	<b>Straight Connector 22/10 Polysulfone MPU</b> For: AH,AK,AL Autoclavable at 134°C Material: Polysulfone 10 mm male to 22 mm female	
156209	<b>Straight Connector 15M/22F MPU</b> For: AI,AH,AK,AL,AN,+nCPAP Autoclavable at 134°C	

## Infant Flow™ nCPAP Generator

Part no.	Description	
7772020LP	Infant Flow™ LP Generator Kit, including generator, sizing guide, and small/medium/large nasal prongs. Single patient use. SPU 10/box	
7772000LP	Infant Flow™ LP Generator Kit, including generator, sizing guide, no nasal prongs. Single patient use. SPU 10/box	
7772011	Infant Flow™ LP nCPAP Generator Circuit Kit (F&P MR 730 & 850 Humidifiers), 10 / case	
7772021	Infant Flow™ LP nCPAP Generator Circuit Kit (F&P MR 730 & 850 Humidifiers), 20 / case	

## Infant Flow™ LP Headgear

Part no.	Description
777040XS	Infant Flow™ LP Headgear, Extra Small, 17-21 cm, 10/package
777040S	Infant Flow™ LP Headgear, Small, 21-26 cm, 10/package
777040SM	Infant Flow™ LP Headgear, Small/Med, 24-28 cm, 10/package
777040M	Infant Flow™ LP Headgear, Medium, 26-32 cm, 10/package
777040L	Infant Flow™ LP Headgear, Large, 32-37 cm, 10/package
777040XL	Infant Flow™ LP Headgear, Extra Large, 37-42 cm, 10/package



## Infant Flow™ LP Bonnets

Part no.	Description
777010	Infant Flow™ LP Bonnet, Size 000, 18-20 cm, 10/package
777012	Infant Flow™ LP Bonnet, Size 00, 20-22 cm, 10/package
777014	Infant Flow™ LP Bonnet, Size 0, 22-24 cm, 10/package
777016	Infant Flow™ LP Bonnet, Size 1, 24-26 cm, 10/package
777018	Infant Flow™ LP Bonnet, Size 2, 26-28 cm, 10/package
777020	Infant Flow™ LP Bonnet, Size 3, 28-30 cm, 10/package
777022	Infant Flow™ LP Bonnet, Size 4, 30-32 cm, 10/package
777024	Infant Flow™ LP Bonnet, Size 5, 32-34 cm, 10/package
777026	Infant Flow™ LP Bonnet, Size 6, 34-36 cm, 10/package
777028	Infant Flow™ LP Bonnet, Size 7, 36-38 cm, 10/package



## Infant Flow™ LP Nasal Prongs

Part no.	Description
777000XS	Nasal Prings, Extra Small, 10 / box
777000S	Nasal Prings, Small, 10 / box
777000M	Nasal Prings, Medium, 10 / box
777000L	Nasal Prings, Large, 10 / box
777000XL	Nasal Prings, Extra Large, 10 / box



## Infant Flow™ LP Nasal Mask

Part no.	Description
777002XS	Nasal Mask, Extra Small, 10 / box
777002S	Nasal Mask, Small, 10 / box
777002M	Nasal Mask, Medium, 10 / box
777002L	Nasal Mask, Large, 10 / box
777002XL	Nasal Mask, Extra Large, 10 / box







## High Flow Nasal Cannula



Part no.	Description
155505	Nasal Cannula O <sub>2</sub> Therapy XXS SPU 50/Box For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy
155515	Nasal Cannula O <sub>2</sub> Therapy XS SPU 50/Box For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy
155525	Nasal Cannula O <sub>2</sub> Therapy S SPU 50/Box For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy
155535	Nasal Cannula O <sub>2</sub> Therapy M SPU 50/Box For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy
155545	Nasal Cannula O <sub>2</sub> Therapy L SPU 50/Box For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy





## SpO<sub>2</sub> Option Masimo

Part no.	Description	
152403	<p><b>LNCS NEO PT SpO<sub>2</sub> Probe, 18" SPU, 20/Box</b>  <b>For: AI,AH,AK,AL,AN,AT</b></p> <p>Patient weight &lt;1 kg                      Single Patient Use                      Latex free                      Masimo provides an adapter from LNCS to RD sensors</p>	
152404	<p><b>LNCS NEO SpO<sub>2</sub> Probe, 18" SPU, 20/Box</b>  <b>For: AI,AH,AK,AL,AN,AT</b></p> <p>Patient weight &lt;3 kg or &gt;40 kg                      Single Patient Use                      Latex free</p>	
152405	<p><b>LNCS INF SpO<sub>2</sub> Probe, 18" SPU, 20/Box</b>  <b>For: AI,AH,AK,AL,AN,AT</b></p> <p>Patient weight 3 kg to 20 kg                      Single Patient Use                      Latex free</p>	
152406	<p><b>Masimo SET, SpO<sub>2</sub> LEMO-HPLPM Module, Set</b>  <b>For: AI,AL,AN,AT</b></p> <p>With LEMO connector, cable length: 182 cm                      Including:                      - 1 pc x LEMO-HPLPM Module                      - 1 pc x each LNCS sensor as a sample:                      LNCS NEO PT, LNCS NEO, LNCS INF</p>	




## etCO<sub>2</sub> Option Microstream

Part no.	Description	
152306	<p><b>Filterline H Set, Airway Adapter, 25/Box</b>  <b>For: AI,AH,AK,AL</b></p>	
152307	<p><b>MicroPod LEMO, Microstream etCO<sub>2</sub> Unit</b>  <b>For: AI,AH,AK,AL</b></p> <p>Unit only, no mounting clamp</p>	

## CO<sub>2</sub> Option Microstream (cont.)



Part no.	Description	
152308	<p><b>Capnoline H INF/NEO, Nasal, 25/Box</b>  <b>For: AI,AH,AK,AL</b></p> <p>Aerogen SOLO T-piece 15mm adapter</p>	
152309	<p><b>MicroPod Clamp Mounting</b>  <b>For: AI,AH,AK,AL</b></p> <p>For mounting on pole or standardized rail</p>	

## etCO<sub>2</sub> Option Lo-Flow






Part no.	Description	
152330	<p><b>Lo-Flo Sidestream etCO<sub>2</sub> Set, Philips</b>  <b>For: AI,AH,AK,AL</b></p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x etCO<sub>2</sub> module REF 152331</li> <li>- 1 pc x etCO<sub>2</sub> mount REF 152332</li> </ul>	
152333	<p><b>etCO<sub>2</sub> Airway Adapter Set H, 10/Box</b>  <b>For: AI,AH,AK,AL</b></p>	
152334	<p><b>etCO<sub>2</sub> Nasal Cannula - Infant, 10/Box</b>  <b>For: AI,AH,AK,AL</b></p>	
152335	<p><b>etCO<sub>2</sub> Nasal Cannula - Paediatric, 10/Box</b>  <b>For: AI,AH,AK,AL</b></p>	



## Aerogen Solo USB Version

Part no.	Description	
152220	<b>Aerogen USB Nebulizer Set</b> <b>For: AI,AH,AK,AN</b> Consisting of: <ul style="list-style-type: none"><li>- 1 pc x Aerogen Solo USB Controller REF 152221</li><li>- 2 pc x Aerogen Solo Nebuliser REF 152222</li><li>- 1 pc x Aerogen Solo T-piece 15mm REF 152223</li><li>- 1 pc x Aerogen Solo T-piece 22mm REF 152224</li></ul>	
152222	<b>Aerogen Solo Nebulizer SPU, 5/Box</b> <b>For: AI,AH,AK,AL,AN</b> Aerogen Solo Nebulizer	
152223	<b>Aerogen Solo T-piece 15mm SPU 10/Box</b> <b>For: AI,AH,AK,AL,AN</b> Aerogen SOLO T-piece 15mm adapter	
152224	<b>Aerogen Solo T-piece 22mm SPU, 10/Box</b> <b>For: AI,AH,AK,AL,AN</b> Aerogen SOLO T-piece 15mm adapter	

## Aerogen PRO Device

Part no.	Description	
152210.AU	<p><b>Aerogen Pro Nebulizer Kit, AU</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For Australia (AU)</p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Nebuliser with filler cab</li> <li>- 1 pc x T-piece with plug</li> <li>- 1 pc x Aerogen Pro controller</li> <li>- 1 pc x Controller cable</li> <li>- 1 pc x AC/DC adapter</li> <li>- 1 pc x Universal mounting bracket &amp; equipment mountadapter</li> </ul>	
152210.EU	<p><b>Aerogen Pro Nebulizer Kit, EU</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For European Union (EU)</p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Nebuliser with filler cab</li> <li>- 1 pc x T-piece with plug</li> <li>- 1 pc x Aerogen Pro controller</li> <li>- 1 pc x Controller cable</li> <li>- 1 pc x AC/DC adapter</li> <li>- 1 pc x Universal mounting bracket &amp; equipment mountadapter</li> </ul>	
152210.UK	<p><b>Aerogen Pro Nebulizer Kit, UK</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For United Kingdom (UK)</p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Nebuliser with filler cab</li> <li>- 1 pc x T-piece with plug</li> <li>- 1 pc x Aerogen Pro controller</li> <li>- 1 pc x Controller cable</li> <li>- 1 pc x AC/DC adapter</li> <li>- 1 pc x Universal mounting bracket &amp; equipment mountadapter</li> </ul>	
152211	<p><b>Aerogen Pro Nebulizerfiller with Cap MPU</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Autoclavable</p>	
152212	<p><b>Aerogen Pro Adapter, NEO/PED, MPU 5/Box</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Paediatric T-adapter &amp; neonate adapter kit</p>	



## Aerogen PRO Device (cont.)

152213

**Aerogen Pro Filler Cap Silic, MPU, 5/Box**  
**For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy**







Aerogen Pro Ped/Neo nebulizer & T-piece kit  
Autoclavable







## Carts

Part no.	Description	
401100	<p><b>Cart for fabian™ HFO</b> <b>For: AI, AH, AL, AK</b></p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Mounting plate for HFO REF 401101</li> <li>- 1 pc x Cart HFO REF 400001</li> </ul> <p>Note: Mounting plate and cart cannot be ordered separately Replaces REF 401210 when stock is empty</p>	
401200	<p><b>Cart for fabian™ EVO Cart</b> <b>For: AN,AT</b></p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Mounting plate for EVO REF 401202</li> <li>- 1 pc x Cart EVO REF 400002</li> </ul> <p>Note: Mounting plate and cart cannot be ordered separately Replaces REF 401210 when stock is empty</p>	



## Air Compressors

Part no.	Description	
490052.EU	<p><b>Compressor Silent 230V/50Hz</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Device on wheels                      Including:                      - 1 pc x 2 m tube                      - 1 pc x Power cord (REF 497180.XX)                      Please choose a power cord when you place the order!</p>	
490058	<p><b>Compressor Silent 230V/60Hz</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For use with Smart Trolley REF 493000                      Including:                      - 1 pc x 2 m tube                      - 1 pc x Power cord (REF 497180.XX)                      Please choose a power cord when you place the order!</p>	
490059.EU	<p><b>Compressor Silent 230V/50Hz</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For use with Smart Trolley REF 493000                      Including:                      - 1 pc x 2 m tube                      - 1 pc x Power cord (REF 497180.XX)                      Please choose a power cord when you place the order!</p>	
490069.US	<p><b>Compressor Silent 115V/60Hz</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For use with Smart Trolley REF 493000                      Including:                      - 1 pc x 2 m tube                      - 1 pc x Power cord (REF 497180.XX)                      Please choose a power cord when you place the order!</p>	
491001.EU	<p><b>Trolley &amp; Compressor 230/50 Set HFO</b>  <b>For: AI,AH,AK,AL</b></p> <p>Consisting of:                      - 1 pc x Smart trolley for compressor mounting REF 493000                      - 1 pc x Mounting plate only for HFO REF 401101                      - 1 pc x Compressor silent 230V/50Hz REF 490059.EU</p>	
491001.SA	<p><b>Trolley &amp; Compressor 230/60 Set HFO</b>  <b>For: AI,AH,AK,AL</b></p> <p>Consisting of:                      - 1 pc x Smart trolley for compressor mounting REF 493000                      - 1 pc x Mounting plate only for HFO REF 401101                      - 1 pc x Compressor silent 230V/60Hz REF 490058</p>	



## Air Compressors (cont.)

Part no.	Description	
491001.US	<p><b>Trolley &amp; Compressor 115/60 Set HFO</b>  <b>For: AI,AH,AK,AL</b></p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Smart trolley for compressor mounting REF 493000</li> <li>- 1 pc x Mounting plate only for HFO REF 401101</li> <li>- 1 pc x Compressor silent 115V/60Hz REF 490069.US</li> </ul>	
491002.EU	<p><b>Trolley &amp; Compressor 230/50 Set EVO</b>  <b>For: AN,AT</b></p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Smart trolley for compressor mounting REF 493000</li> <li>- 1 pc x Mounting plate only for evolution REF 401202</li> <li>- 1 pc x Compressor silent 230V/50Hz REF 490059.EU</li> </ul>	
491002.SA	<p><b>Trolley &amp; Compressor 230/60 Set EVO</b>  <b>For: AN,AT</b></p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Smart trolley for compressor mounting REF 493000</li> <li>- 1 pc x Mounting plate only for evolution REF 401202</li> <li>- 1 pc x Compressor silent 230V/60Hz REF 490058</li> </ul>	
491002.US	<p><b>Trolley &amp; Compressor 115/60 Set evo</b>  <b>For: AN,AT</b></p> <p>Consisting of:</p> <ul style="list-style-type: none"> <li>- 1 pc x Smart Trolley for compressor mounting REF 493000</li> <li>- 1 pc x Mounting plate only for evolution REF 401202</li> <li>- 1 pc x Compressor silent 115V/60Hz REF 490069.US</li> </ul>	

## Carts Accessories

Part no.	Description	
400011.01	<p><b>IV-pole with DIN Clamp</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS3</b></p> <p>Short IV-pole, length 1.2 m                      With height adjustment                      Universal clamp for standard bar (25 x 10) mm and round pole 25 mm</p>	
410001	<p><b>Circuit Support Arm with DIN Rail Clamp</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS3</b></p> <p>One knob for fixation and adjustment                      DIN rail clamp bracket</p>	
400013	<p><b>Basket, AMS Cart &amp; Light Cart</b>  <b>for Tolley REF 401100 or 401200</b></p>	
400014	<p><b>Bottle Holder, AMS Cart &amp; Light Cart</b>  <b>for Tolley REF 401100 OR 401200</b></p>	




## Carts Spares

Part no.	Description	
401101	<p><b>Mounting Plate only for HFO Cart</b>  <b>For: AI,AH,AK,AL</b></p> <p>Material: Stainless steel, powder-coated                      Inclusive:                      - 2 pcs x DIN rails (square tube 25/10)                      - Fixing material</p>	
401202	<p><b>Mounting Plate only for EVO Cart</b>  <b>For: AN,AT,+nCPAP,Therapy</b></p> <p>Material: Stainless steel, powder-coated                      Inclusive:                      - 2 pcs x DIN rails (square tube 25/10)                      - Fixing material</p>	






## Extended Warranty Program

Part no.	Description
111191	<b>fabian™ HFO, Warranty Extension +1 Year</b> For: AI,AH,AK,AL
111192	<b>fabian™ HFO, Warranty Extension +2 Years</b> For: AI,AH,AK,AL
111193	<b>fabian™ HFO, Warranty Extension +3 Years</b> For: AI,AH,AK,AL
121191	<b>fabian™ Therapy evolution Warranty +1 Year</b> For: AT
121192	<b>fabian™ Therapy evolution Warranty +2 Years</b> For: AT
121193	<b>fabian™ Therapy evolution Warranty +3 Years</b> For: AT
122191	<b>fabian™ +nCPAP evolution Warranty +1 Year</b> For: AN
122192	<b>fabian™ +nCPAP evolution Warranty +2 Years</b> For: AN
122193	<b>fabian™ +nCPAP evolution Warranty +3 Years</b> For: AN


## Air and O<sub>2</sub> Sets

Part no.	Description	
499090	<b>Filterset for Medical Grade Air-NIST</b> For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS To be mounted on the trolley or on the wall, screws are included	
499091.EU	<b>Air Water Trap NIST 90°</b> For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS3 Water trap with female NIST adapter on ventilator side and male NIST adapter on input side	
499200.DE	<b>Air/O<sub>2</sub> Hose Set 5M BL/WH, DIN120 - NIST</b> For: MS3 Consisting of: - 1 pc x AIR hose, BK-WH, 5m, NIST-DIN - 1 pc x O <sub>2</sub> hose, WH, 5m, NIST-DIN	









## Air and O<sub>2</sub> Sets (cont.)

Part no.	Description	
499230.BS	<p><b>Air/O<sub>2</sub> Hose Set, 3m, BL/WH, NIST-BS5682</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Consisting of:                      - 1 pc x Air hose, BK-WH, 3 m, NIST-BS5682                      - 1 pc x O<sub>2</sub> hose, WH, 3 m, NIST-BS5682</p>	
499230.DE	<p><b>Air/O<sub>2</sub> Hose Set, 3m, NIST-DIN</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Consisting of:                      - 1 pc x Air hose, BK-WH, 3 m, NIST-DIN                      - 1 pc x O<sub>2</sub> hose, WH, 3 m, NIST-DIN</p>	
499301.US	<p><b>Air Water Trap NIST 90° Male #A3/DISS</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>For USA, Canada and Japan</p>	
499391.US	<p><b>Air Adapter, NIST-f / DISS-m, 90°</b>  <b>For: AI,AH,AK,AL,AN,AT,MS3</b></p>	
499392.US	<p><b>O<sub>2</sub> Adapter, NIST-f / DISS-m, 90°</b>  <b>For: AI,AH,AK,AL,AN,AT,MS3</b></p> <p>For USA, Canada and Japan</p>	

## Powercords Angled

Part no.	Description	
497090.AU	<p><b>Powercord AUS, 90° (C13 / AS3112), 2.5 m</b>  <b>For: AI,AH,AK,AN,AT,+nCPAP,Therapy,MS3</b></p> <p>Australia                      Cannot be ordered separately!                      For Compressor, Humidifier, InfantView</p>	
497090.CH	<p><b>Powercord CH, 90° (C13 / T12), 2.5 m</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS3</b></p> <p>Switzerland                      Cannot be ordered separately!</p>	
497090.EU	<p><b>Powercord EU, 90° (CEE 7/7), 2.5 m</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS3</b></p> <p>SCHUKO / Europe                      Cannot be ordered separately!</p>	

## Powercords Angled (cont.)

Part no.	Description	
497090.JP	<p><b>Powercord JP, 90° (C13 / T12), 2 m</b>  <b>AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS3</b></p> <p>USA/Japan            Cannot be ordered separately!</p>	
497090.UK	<p><b>Powercord UK, 90° (C13 / BS1363), 2.5 m</b>  <b>AI,AH,AK,AL,AN,AT,+nCPAP,Therapy,MS3</b></p> <p>United Kingdom            Cannot be ordered separately!</p>	
497180.01US	<p><b>Powercord US JP, Type-B - NEMA</b>  <b>For: AI,AH,AK,AL,AN,AT</b></p> <p>USA/Japan            Cannot be ordered separately!</p>	
497180.02US	<p><b>Powercord US, Hosp. Grade 10A, Trnsp 3 m</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>USA            Cannot be ordered separately!</p>	
497180.AU	<p><b>Powercord AUS, (C13 / AS3112), 2.5 m</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Australia            Cannot be ordered separately!</p>	
497180.CH	<p><b>Powercord CH, (C13 / T12), 2.5 m</b>  <b>For: AI,AH,AK,AL,AN,AT</b></p> <p>Switzerland            Cannot be ordered separately!</p>	
497180.EU	<p><b>Powercord EU (CEE 7/7), 2.5 m</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>SCHUKO / Europe            Cannot be ordered separately!</p>	
497180.UK	<p><b>Powercord UK, (C13 / BS1363), 2.5 m</b>  <b>For: AI,AH,AK,AL,AN,AT,+nCPAP,Therapy</b></p> <p>Straight, United Kingdom            Cannot be ordered separately!</p>	



## GLOBAL HEADQUARTERS

Vyaire Medical, Inc.  
26125 North Riverwoods Blvd  
Mettawa, IL 60045  
USA

[vyaire.com](http://vyaire.com)

### For International distribution.

© 2021 Vyaire. Vyaire, the Vyaire logo and all other trademarks or registered trademarks are property of Vyaire Medical, Inc., or one of its affiliates. VYR-GBL-2100261 v2.0

