

Product Identification	
Product family	Sterile suture sets
Reference	10608088
Name	Sutureringsæt
Sterilization	Ethylene oxide
Brand	Sterisets®



Material			
Code	Designation	Qty	Composition
213004INT	Tray PP 180 x 95 x 25 mm 2 compartments transp.	1	Polypropylene
56405PC5	Swab NW 4 ply 30 gr 10 x 10 cm counted per 5	3	Non woven 30 g/m2
56201	Ball gauze 17 thr. 12 x 15 cm Ø 30 mm with elastic ring	5	Absorbent gauze 100% cotton
570743	Needleholder Crile-Murray 15cm spec, design fine tip no card Blue	1	AISI-410 Stainless steel
570034	Forceps dressing delicate Adson 1 x 2 teeth smooth 12 cm w/c	1	AISI-410 Stainless steel
570059	Scissors iris delicate S/S CVD 11 cm w/c	1	AISI-420 Stainless steel
1003INT	Wrapping field PE/Tissue 50 x 50 cm green	1	Polyethylene/tissue 10 µm

Packaging	
Type	Sterile single pack
Package weight	114g
Package dimensions	260 x 130 x 20mm
Quantity per inner box	28pcs
Quantity transport carton	112pcs

### Hazards Identification

Latex	No
Medicinal substances	No
Animal tissue	No
Toxic or dangerous materials	No

### Shelf life

5 years from the date of manufacture, unless any component has a shorter shelf-life.

### Storage conditions

Keep away from sunlight and keep dry.

### Environmental and disposal considerations

Not regarded as dangerous for the environment. Dispose of material according to local, regional and national regulations.

### Compliance Regulations & Standards

Regulation (EU) 2017/745 (as applicable), Directive 93/42/EEC, ISO 13485, ISO 14971, ISO 15223-1, ISO 10993-1, ISO 11737, ISO 11135.

