

Safety Scalpel

Stainless Steel blade


Features

- It is a safe alternative to the traditional scalpel and fulfills the requirement to reduce sharp injuries.
- Ribbed on both side of the sheath which provide extra gripping for left handed as well as right handed surgeon.
- One-handed activation allowing the surgeon to focus on the surgical site.
- Audible clicks makes the scalpel more secure and protective.
- "Used & Un-Used" indications prevent contamination and transmission of blood borne disease.
- Permanent safety lock prevents accidental cuts.
- Latex Free.
- Measurement scale 0 to 5 cm is embossed on the handle.
- Ergonomically designed safety scalpel.
- ABS plastic handle and sheath gives extra strength to the product.
- Universal size handle that is ideal for kits, trays & packs.


Details:

Product Name:	Safety Scalpel
Size of Inner Box (cm):	10 Pcs
Unit/Box:	21.5 x 7.8 x 3.5
Outer Boxes:	5 Boxes x 10 Units
Export carton:	5 Boxes x 20 Outer
Export Carton size (cm):	48 x 46.5 x 38.5
Per 20 feet container:	300 Carton


Size	Product Code	Fitment No	Size	Product Code	Fitment No
10	PSF10	3	21	PSF21	4
10A	PSF10A	3	22	PSF22	4
11	PSF11	3	23	PSF 23	4
11K	PSF11K	3	24	PSF24	4
13	PSF13	3	25	PSF25	4
15	PSF15	3	36	PSF36	4
15B	PSF15B	3	MSC	PSFMSC	3
15C	PSF15C	3			
17	PSF17	4			
18	PSF18	4			
19	PSF19	4			
20	PSF20	4			




SIZE 10



SIZE 11




SIZE 12



SIZE 15




SIZE 20



SIZE 21



SIZE 22



SIZE 23

