

Feeding Tubes

Nutricare Flo is an open-ended polyurethane tube with clearly numbered cm markings. It is fitted with an ENFit® (ISO 80369-3) connector

Enteral Nutricare Flo

- Made of PU
- DEHP & latex free
- Male ENFit® connector
- Radiopaque stripe
- Clearly numbered every cm
- Polyurethane feeding tube with a rounded open distal end
- *“Open ended feeding tubes have been reported to lead to a reduction in bacterial colonisation”¹*

French Size	Length	Priming Volume
4Fr	50cm	0.25ml
	75cm	0.75ml
5Fr	50cm	0.5ml
	75cm	0.75ml
6Fr	50cm	0.70ml
	75cm	1.1ml
8Fr	50cm	1.25ml
	75cm	1.85ml

		25 pcs / box Product Code:	NHSSC:
4Fr	40cm	NCF04.40	FWM4381
	50cm	NCF04.50	TBC
	75cm	NCF04.75	TBC
5Fr	40cm	NCF05.40	FWM4382
	50cm	NCF05.50	TBC
	60cm	NCF05.60	FWM4373
	75cm	NCF05.75	TBC
	90cm	NCF05.90	FWM4379
6Fr	125cm	NCF05.125	TBC
	40cm	NCF06.40	FWM4380
	50cm	NCF06.50	TBC
	60cm	NCF06.60	FWM4374
	75cm	NCF06.75	TBC
	90cm	NCF06.90	FWM4383
8Fr	125cm	NCF06.125	TBC
	40cm	NCF08.40	FWM4376
	50cm	NCF08.50	TBC
	60cm	NCF08.60	FWM4372
	75cm	NCF08.75	TBC
10Fr	90cm	NCF08.90	FWM4378
	125cm	NCF08.125	FWM4377
	75cm	NCF10.75	TBC
12Fr	110cm	NCF10.110	TBC
	110cm	NCF12.110	TBC

¹ Enteral Feeding Tube Design and Differential Bacterial Overgrowth: An In Vitro Comparison. Amy Presti, MD, MSCR; Ruth Snyder, NNP, BC: Neonatal INTENSIVE CARE Vol. 26 No. 4 n July-August 2013

