

# Needles & Syringes

For an extra secure connection the threaded tip provides a locking fit. Compatible with a variety of needles, cannulas and other IV devices

## Luer Lock Syringe

### Features & Benefits

- Bold easy to read graduations
- Available in a full range of sizes
- Smooth plunger motion
- Latex & PVC free
- Integrated plunger backstop
- Highly transparent barrel

**60ml**  
30 pcs / box  
Product Code: S0-600



**30ml**  
50 pcs / box  
Product Code: S0-300



**20ml**  
100 pcs / box  
Product Code: S0-200



**10ml**  
100 pcs / box  
Product Code: S0-100



**5ml**  
100 pcs / box  
Product Code: S0-050




**3ml**  
100 pcs / box  
Product Code: S0-030



### Low Dead Space Spigot Tip Design

**1ml**  
100 pcs / box  
Product Code: S0-010




# Needles & Syringes

A friction fit connection that requires the clinician to insert the tip of the syringe into the needle hub or other attaching device in a push and twist manner. This twist will ensure a secure connection that is less likely to detach

## Luer Slip Syringe

### 1ml Tuberculin Syringe

Historically used to deliver the TB antigen test, Tuberculin syringes are typically 1ml size, and used to deliver very small amounts of medication accurately, through a needle. The GBUK Tuberculin syringe features a dead space reducing spigot in the tip, for extra accuracy



**1ml**  
100 pcs / box  
Product Code: S0-011T

### Low Dead Space Spigot Tip Design

**1ml**  
100 pcs / box  
Product Code: 3PLS001



**3ml**  
100 pcs / box  
Product Code: 3PLS003




**5ml**  
100 pcs / box  
Product Code: 3PLS005




**10ml**  
100 pcs / box  
Product Code: 3PLS010



**20ml**  
50 pcs / box  
Product Code: 3PLS020



**30ml**  
40 pcs / box  
Product Code: 3PLS030



**60ml**  
30 pcs / box  
Product Code: 3PLS060



### Features & Benefits

- Bold easy to read graduations
- Available in a full range of sizes
- Smooth plunger motion
- Latex & PVC free
- Integrated plunger backstop
- Highly transparent barrel

