



# NRFit<sup>®</sup> Syringe

Intended for the administration of medication or anesthetic via the Neuraxial route. ISO 80369-6:2016 is the international standard that specifies requirements for small-bore connectors intended to be used for connections in neuraxial applications

# Lock Syringe



Not made with natural rubber latex



Connectors compliant with ISO 80369-6 are specific to neuraxial devices



NRFit<sup>®</sup> devices from GBUK Healthcare are exclusively branded DASH 6<sup>™</sup>



The term "NRFit<sup>®</sup>" refers to neuraxial devices that are ISO 80369-6 compliant

Lock Syringe



## Neuraxial connectors are to be branded NRFit

Where practicable these devices will for easy identification be coloured yellow

<b>1ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15011010</b>
BNS (1000pcs / box)	<b>15010010</b>






<b>3ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15011030</b>
BNS (500pcs / box)	<b>15010030</b>

<b>5ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15010050</b>
BNS (500pcs / box)	<b>15011050</b>

<b>10ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15011100</b>
BNS (500pcs / box)	<b>15010100</b>

<b>20ml</b>	Product Code:
Sterile (50 pcs / box)	<b>15011200</b>
BNS (400pcs / box)	<b>15010200</b>

<b>60ml</b>	Product Code:
Sterile (25 pcs / box)	<b>15011600</b>
BNS (200pcs / box)	<b>15010600</b>



# NRFit<sup>®</sup> Lock Syringe





# NRFit<sup>®</sup> Syringe

Intended for the administration of medication or anesthetic via the Neuraxial route. ISO 80369-6:2016 is the international standard that specifies requirements for small-bore connectors intended to be used for connections in neuraxial applications

# Slip Syringe



Not made with natural rubber latex



Connectors compliant with ISO 80369-6 are specific to neuraxial devices



NRFit<sup>®</sup> devices from GBUK Healthcare are exclusively branded DASH 6<sup>™</sup>



The term “NRFit<sup>®</sup>” refers to neuraxial devices that are ISO 80369-6 compliant



## Neuraxial connectors are to be branded NRFit

Where practicable these devices will for easy identification be coloured yellow

<b>1ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15021010</b>
BNS (1000pcs / box)	<b>15020010</b>

<b>3ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15021030</b>
BNS (500pcs / box)	<b>15020030</b>

<b>5ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15021050</b>
BNS (500pcs / box)	<b>15020050</b>

<b>10ml</b>	Product Code:
Sterile (100 pcs / box)	<b>15021100</b>
BNS (500pcs / box)	<b>15020100</b>

<b>20ml</b>	Product Code:
Sterile (50 pcs / box)	<b>15021200</b>
BNS (400pcs / box)	<b>15020200</b>



# NRFit<sup>®</sup> Slip Syringe






# NRFit<sup>®</sup> Accessories


NRFit<sup>®</sup> Accessories have been developed to improve the use of medical devices intended to administer medications to neuraxial sites, wound infiltration anaesthesia delivery, and other regional anaesthesia procedures or to monitor or remove cerebrospinal fluid for therapeutic or diagnostic purposes

# NRFit<sup>®</sup> Accessories

## Blunt Drawing Up Needles

## Extension Sets

 Not made with natural rubber latex

 ISO 80369-6 NRFit<sup>®</sup>



**Blunt Drawing Up Needle 5µm x 50mm**

	Product Code:
Sterile (100 pcs / box)	15291050
BNS (1000pcs / box)	15490050



**Blunt Drawing Up Needle 5µm x 100mm**

	Product Code:
Sterile (100 pcs / box)	15281100
BNS (1000pcs / box)	15240001



**Blunt Drawing Up Needle 5µm x 50mm**

	Product Code:
Sterile (100 pcs / box)	15491050
BNS (1000pcs / box)	



**21G Butterfly Needle One Way Valve**

	Product Code:
Sterile (100 pcs / box)	15361221



**25G Butterfly Needle**

	Product Code:
Sterile (100 pcs / box)	15351225



**27G Butterfly Needle**

	Product Code:
Sterile (100 pcs / box)	15351227




**Extension Set With Clamp 10cm**

	Product Code:
Sterile (25 pcs / box)	15461010
BNS (400pcs / box)	15460010



**Extension Set With Clamp 20cm**

	Product Code:
Sterile (25 pcs / box)	15461020
BNS (400pcs / box)	15460020



**Extension Set With Clamp 125cm**

	Product Code:
Sterile (25 pcs / box)	15461125
BNS (200pcs / box)	15460125

# NRFit<sup>®</sup> Accessories



Scan for more information & order codes

# NRFit® Accessories



NRFit® Accessories have been developed to improve the use of medical devices intended to administer medications to neuraxial sites, wound infiltration anaesthesia delivery, and other regional anaesthesia procedures or to monitor or remove cerebrospinal fluid for therapeutic or diagnostic purposes


# NRFit® Accessories



Not made with natural rubber latex



ISO 80369-6  
NRFit®



Cap & Adaptors



Filters



Stopcock



Drawing Up



Introducer



## Needle Hub Cap

	Product Code:
Sterile (100 pcs / box)	15111001
BNS (1000pcs / box)	15110001



## Syringe Cap

	Product Code:
Sterile (100 pcs / box)	15111001
BNS (1000pcs / box)	15110001



## Syringe to Syringe Adaptor

	Product Code:
Sterile (100 pcs / box)	15171001



## 0.22µm Disc Filter

	Product Code:
Sterile (100 pcs / box)	15231022



## 3-Way Stopcock

	Product Code:
Sterile (100 pcs / box)	15241001
BNS (1000pcs / box)	15240001




## Epidural Flat Filter

	Product Code:
Sterile (100 pcs / box)	15311022
BNS (250pcs / box)	15310022



## Epidural Catheter Feeder

	Product Code:
BNS (1000pcs / box)	15200001



## Tuohy Borst Adaptor

	Product Code:
BNS (1000pcs / box)	15300001



## Drawing Up Filter Straw, 5µm, 100mm

	Product Code:
Sterile (100 pcs / box)	15281100
BNS (1000pcs / box)	15280100



## Introducer Needle 18G x 38mm

	Product Code:
Sterile (400 pcs / box)	15090002

# NRFit® Accessories

