

Needles & Syringes

Hypodermic Non-Safety Needle



Wide Range of Gauges & Lengths

each needle hub is colour coded for easy gauge identification






Extremely Smooth Surface

with light silicone coating, for minimal pain on insertion

Additional full range of sizes also available (18G-31G/ 0.5-1/5"). Please ask your Area Business Manager for further details

Needle with thin wall technology. 18G-25G. Reduces patient discomfort and risk of haemolysis¹


1. Amy Mouser, et al, "Evaluation of Pain and Specimen Quality by Use of a Novel 25-Gauge Blood Collection Set With Ultra-Thin Wall Cannula and 5-Bevel Tip Design," *The Journal of Applied Laboratory Medicine*, 2017, 2(2):201-210

18G	100 pcs / box Product Code:	
40mm (1.5 inch)	HN181500	
50mm (2 inch)	HN182000	
19G	100 pcs / box Product Code:	
30mm (1.25 inch)	HN191250	
38mm (1.5 inch)	HN191500	
40mm (1.5 inch)	HN191500	
50mm (2 inch)	HN192000	
20G	100 pcs / box Product Code:	
25mm (1 inch)	HN201000	
38mm (1.5 inch)	HN201500	
50mm (2 inch)	HN202000	
21G	100 pcs / box Product Code:	
25mm (1 inch)	HN211000	
38mm (1.5 inch)	HN211500	
50mm (2 inch)	HN212000	
22G	100 pcs / box Product Code:	
30mm (1.25 inch)	HN221250	
50mm (2 inch)	HN222000	
23G	100 pcs / box Product Code:	
25mm (1 inch)	HN231000	
32mm (1.25 inch)	HN231250	
25G	100 pcs / box Product Code:	
16mm (0.625 inch)	HN250625	
25mm (1 inch)	HN251000	
26G	100 pcs / box Product Code:	
10mm (0.375 inch)	HN260375	
13mm (0.5 inch)	HN260500	
25mm (1 inch)	HN261000	
27G	100 pcs / box Product Code:	
13mm (0.5 inch)	HN270500	
19mm (0.75 inch)	HN270750	
38mm (1.5 inch)	HN271500	
30G	100 pcs / box Product Code:	
13mm (0.5 inch)	HN300500	



Hypodermic Safety Needle

Features



Wide Range of Gauges and Lengths

each needle hub is colour coded for easy gauge identification



Extremely Smooth Surface

with light silicone coating for minimal pain on insertion



Safety Needle

safety needle with one handed intuitive activation and audible click

18G	100 pcs / box Product Code:
38mm (1.5 inch)	SHN181500
50mm (2 inch)	SHN182000

19G	100 pcs / box Product Code:
25mm (1 inch)	SHN191000
38mm (1.5 inch)	SHN191500
50mm (2 inch)	SHN192000

20G	100 pcs / box Product Code:
25mm (1 inch)	SHN201000
38mm (1.5 inch)	SHN201500
50mm (2 inch)	SHN202000



Safety Cover

Completely houses the entire needle for total security

21G	100 pcs / box Product Code:
25mm (1 inch)	SHN211000
38mm (1.5 inch)	SHN211500
50mm (2 inch)	SHN212000

22G	100 pcs / box Product Code:
25mm (1 inch)	SHN221000
38mm (1.5 inch)	SHN221500
50mm (2 inch)	SHN222000

23G	100 pcs / box Product Code:
25mm (1 inch)	SHN231000
32mm (1.25 inch)	SHN231250
40mm (1.5 inch)	SHN231500

25G	100 pcs / box Product Code:
16mm (0.625 inch)	SHN250625
25mm (1 inch)	SHN251000
38mm (1.5 inch)	SHN251500

26G	100 pcs / box Product Code:
13mm (0.5 inch)	SHN260500
16mm (0.625 inch)	SHN260625

27G	100 pcs / box Product Code:
16mm (0.625 inch)	SHN270625
25mm (0.75 inch)	SHN270750
32mm (1.25 inch)	SHN271250

30G	100 pcs / box Product Code:
13mm (0.5 inch)	SHN300500

