








IV3000 range

Dressing	Size (cm)	S&N code	NHS code	Quantity	Recommended indications
	5 x 6	66004011	ELW011	100	Paediatric/Ported
	6 x 7	4007	ELW032	100	Peripheral
	6 x 7	59410082	ELW105	100	Peripheral
	7 x 9	4006	ELW046	100	Ported/Peripheral
	9 x 12	66004009	ELW112	50	Central/Jugular
	10 x 12	4008	ELW054	50	Central
	10 x 12	59410882	ELW106	50	Central
	11 x 14	66800512	ELW363	25	PICC/Midline
	10 x 20	4649	ELW099	50	Central/Epidural
OPSITE^o Post-Op					
	5 x 6.5	66000708	ELW052	100	Post removal of IV line
ACTICOAT^o Site					
	2.5	66001788	ELY281	10	Central/Jugular/PICC

Dressing		Size's available/code	Carton size
IV3000	Ported 	5cm x 6cm (66004011) 9cm x 12cm (66004009) 7cm x 9cm (4006) 11cm x 14cm (66800512)	100 50 100 25
	Non-ported 	6cm x 7cm (4007) 10cm x 12cm (4008)	100 50
IV3000 Frame Delivery		6cm x 7cm (59410082) 10cm x 12cm (59410882)	100 50
IV3000 Standard		6cm x 8cm (4923) 6cm x 8.5cm (4924) 10cm x 14cm (4925) 10cm x 20cm (4649) 10cm x 14cm (4973)	100 100 10 50 50
IV3000 Infusion Set		6cm x 7cm ported (66800786)	30

IV3000°

Moisture responsive catheter dressing

Indications

For dressing central line, peripheral, PICC, epidural and paediatric catheters.

Contra-indications

The dressings should be used on IV or epidural sites. If a film dressing for open wounds, such as pressure sores, burns and donor sites is required, OPSITE® FLEXIGRID® should be used.

IV3000 I-Hand only:

The securing strips and documentation label are not indicated as a replacement for sutures or other primary wound closure methods.



Aquillant Medical

IV3000 is distributed in the UK by Aquillant Medical

S&N Code	Description	Pack Content	NHSSC Code	PIP Code
66004011*	5cm x 6cm ported paediatric catheter	100	ELW011	-
4006*	7cm x 9cm ported peripheral catheters	100	ELW046	276-0783
66004009*	9cm x 12cm central/jugular catheters	50	ELW112	320-2827
4007*	6cm x 7cm peripheral catheters	100	ELW032	238-0269
66004005*	10cm x 12cm central line catheters	10	-	305-7411
4008*	10cm x 12cm central line catheters	50	ELW054	238-0285
Frame Delivery Application				
59410082*	6cm x 7cm peripheral catheters	100	ELW105	289-6710
59410882*	10cm x 12cm central catheters	50	ELW106	289-6702
IV3000 Orange Handle application				
4923	6cm x 8cm ported peripheral catheters	100	ELW037	238-0251
4924	6cm x 8.5cm peripheral catheters	100	ELW039	320-2819
4925	10cm x 14cm central line catheters	10	ELW041	238-0277
4649	10cm x 20cm epidural catheters	50	ELW099	-

* 2 sterile securing strips and date label included