

KLINIDERM - Foam Silicone Border

Description

Adhesive foam silicone dressing. For the treatment of chronic and acute wounds. The foam dressing is soft, comfortable and absorbent. Creates the ideal environment for moist wound healing and avoids adherence to the wound surface. The skin-friendly, adhesive dressing enables an optimal fixation. Medical device. Sterile if pouch and seal are intact.

The range offers a wide variety of shaped products for different application areas. The heel dressing can be placed on the heel, whereas the sacrum dressing ensures optimal comfort on the sacral area.









To ensure optimal wear comfort to the patient, we added a new size to our Kliniderm® Foam Silicone range. The Kliniderm® Foam Silicone Multisite can be used on parts of the body where extra conformability is desired. The design is specially designed to stay comfortably on the wound even on the difficult to dress parts of the body.

Product information

In compliance with: MDD (EU) 93/42/EEC Class IIb



Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3016535	KLINIDERM FOAM SILICONE BORDER 7,5X7,5CM STERIEL 40514828	5 pcs.	100 pcs.			40514828
	3016536	KLINIDERM FOAM SILICONE BORDER 10X10CM STERIEL 40514829	5 pcs.	100 pcs.			40514829
	3016537	KLINIDERM FOAM SILICONE BORDER 12,5X12,5CM STERIEL 40514830	5 pcs.	100 pcs.			40514830
	3016538	KLINIDERM FOAM SILICONE BORDER 15X15CM STERIEL 40514831	5 pcs.	100 pcs.			40514831
	3016539	KLINIDERM FOAM SILICONE BORDER 10X20CM STERIEL 40514832	5 pcs.	100 pcs.			40514832
	3016540	KLINIDERM FOAM SILICONE BORDER 10X30CM STERIEL 40514833	5 pcs.	100 pcs.			40514833
	3016541	KLINIDERM FOAM SILICONE BORDER 15X20CM STERIEL 40514834	5 pcs.	100 pcs.			40514834
	3016544	KLINION KLINIDERM FOAM SILICONE SACRUM BORDER WOUND DRESSING 22.50 X 22.50 CM REF 40514837 STERIEL	5 pcs.	100 pcs.			40514837

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3016546	KLINION KLINIDERM FOAM SILICONE SACRUM BORDER WOUND DRESSING 15 X 15 CM REF 40514839	5 pcs.	100 pcs.			40514839
	3037336	KLINION KLINIDERM FOAM SILICONE HEEL BORDER PRE-SHAPED WOUND DRESSING FOR HEEL 20 X 20.8 CM STERILE REF 40514836	5 pcs.	100 pcs.			40514836
	3067525	KLINIDERM FOAM SIL BRD MULTISITE 20X20CM 40514822 STERIEL	5 pcs.	100 pcs.			40514822