



Needle Hub Cap

	Product Code:	NHSSC:
Sterile (100 pcs / box)	15111001	FWF149
BNS (1000 pcs / box)	15110001	TBC



Syringe Cap

	Product Code:	NHSSC:
Sterile (100 pcs / box)	15041001	FWC489
BNS (1000pcs / box)	15040001	TBC



Syringe to Syringe Adaptor

	Product Code:	NHSSC:
Sterile (100 pcs / box)	15171001	FWF150



0.22µm Disc Filter

	Product Code:	NHSSC:
Sterile (100 pcs / box)	15231022	FWC490



3-Way Stopcock

	Product Code:	NHSSC:
Sterile (100 pcs / box)	15241001	FWC491
BNS (1000pcs / box)	15240001	TBC



Epidural Flat Filter

	Product Code:	NHSSC:
Sterile (100 pcs / box)	15311022	FWC493
BNS (250pcs / box)	15310022	TBC



Epidural Catheter Feeder

	Product Code:	NHSSC:
BNS (1000pcs / box)	15200001	TBC



Tuohy Borst Adaptor

	Product Code:	NHSSC:
BNS (1000pcs / box)	15300001	TBC



Drawing Up Filter Straw, 5µm, 100mm

	Product Code:	NHSSC:
Sterile (100 pcs / box)	15281100	FWF151
BNS (1000pcs / box)	15280100	TBC



Introducer Needle 18G x 38mm

	Product Code:	NHSSC:
Sterile (400 pcs / box)	15090002	FTR3228

NRFit Accessories