

Syringes

Enteral ISOSAF Home, reusable enteral syringe with female ENFit® connection (ISO 80369-3). Low dose tip are available

- Silicone O ring
- DEHP free & Latex free
- Clear barrel markings
- Single patient use up to 14 days*
- 60ml with eccentric tip
- Low dose options available in 1, 3 & 5ml syringes
- Full thread

Enteral ISO Saf Home

		100 pcs / box Product Code:	NHSSC:
1ml	Low Dose:	ELDISO-1H	FTR2105
3ml	Low Dose:	ELDISO-3H	FTR2102
5ml	Low Dose:	ELDISO-5H	FTA440
10ml	Standard ENFit Tip:	EISO-10H	FTR2103
20ml	Standard ENFit Tip:	EISO-20H	FTR2104
		50 pcs / box Product Code:	NHSSC:
60ml	Standard ENFit Tip:	EISO-60H	FTR2106

* Depending on syringe contents and washing conditions
- discard if syringe is no longer usable

