

Klinion® nitrile ultra comfort gloves, 100 pcs

Description

The Klinion Personal Protection examination gloves are intended to be worn on the hand for medical purposes and to help prevent contamination between patient and examiner. This information does not reflect the actual duration of protection in the workspace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. This product provides protection against bacteria, fungi and virus. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. The product contains zinc dibutyl dithiocarbamates and should not be used in case of hypersensitivity to these substances. Before usage inspect the gloves for any defect or imperfections. For single use only. Re-use can cause contamination or infection. Dispose according to local procedures and guidelines. Not made with natural rubber latex. Medical device.

Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I & PPE Regulation (EU) 2016/425 Category III






Klinion Nitrile gloves are available in:

- Ultra Comfort (Thin glove in a white or blue color)
- Ultra Safe (Better protection than Ultra Comfort)
- Sensitive (Accelerator-free gloves to reduce the risk on allergies)

The total length of these gloves is minimal 240 mm.



Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3060938	Klinion Personal Protection Examination Glove Nitrile Ultra Comfort Xs Blue 103520	100 pcs.	1000 pcs.			103520
	3060939	Klinion Personal Protection Examination Glove Nitrile Ultra Comfort S Blue 103521	100 pcs.	1000 pcs.			103521
	3060970	Klinion Personal Protection Examination Glove Nitrile Ultra Comfort M Blue 103522	100 pcs.	1000 pcs.			103522
	3060971	Klinion Personal Protection Examination Glove Nitrile Ultra Comfort L Blue 103523	100 pcs.	1000 pcs.			103523
	3060972	Klinion Personal Protection Examination Glove Nitrile Ultra Comfort XI Blue 103524	100 pcs.	1000 pcs.			103524