

CAREFIX® WOUND CARE FIXATION

GENTLE, TIME-SAVING & COST-EFFICIENT SOLUTIONS

TYTEX.COM/WOUNDCAREFIXATION



carefix
by Tytex

CAREFIX® TUBE FIXATION

PRECIT IV & DRESSING FIXATION

Carefix® tubes are designed for secure and comfortable fixation of IV catheters, IV tubes, and dressings on the hand or foot. Their pre-cut, ready-to-use design ensures quick application, making them ideal for emergency rooms and hospital units where time optimization is critical.

By keeping IV catheters securely in place, Carefix tubes help reduce the frequency of dressing changes, lowering infection risks and overall healthcare costs.

YOUR BENEFITS INCLUDE

- **Secure fixation:** Prevents displacement, pulling, or vein irritation of IV catheters and tubes.
- **Reduces dressing changes:** Helps maintain catheter stability, minimizing the need for frequent replacements.
- **Lower infection risk:** Secure fixation reduces exposure and the need for repeated interventions.
- **Cost-effective:** Reduces material waste and overall healthcare costs associated with IV management.
- **Quick & easy application:** Pre-cut design enables fast and hassle-free use, ideal for high-paced hospital environments.



CAREFIX® TUBE FIXATION

TIME-SAVING & COST-EFFICIENT SOLUTION

Selected features

- Secure IV fixation: Prevents displacement of IV catheters, ensuring safety and stability.
- Skin-friendly protection: Eliminates the need for adhesive bandages, protecting sensitive skin.
- Quick & easy application: Pre-cut lengths and practical thumb opening for an optimal fit.
- Flexible & breathable material: Soft, comfortable, and adaptable to different patient needs.
- Customizable fit: Can be trimmed without fraying; color-coded for easy size identification.
- Durable & reusable: Washable, latex-free, and Oeko-Tex certified for patient and caregiver safety.

Composition

70% Polyamide, 21% Elastane, 9% Polyester

Sizes available

S-XL

Packaging

Dispenser boxes



Washing instructions



CAREFIX[®] FINGER / TOE FIXATION

GENTLE & SECURE BANDAGE FIXATION FOR FINGERS AND TOES

The Carefix[®] Finger/Toe Fixation Device is a soft and elastic dressing holder designed to keep surgical dressings and wound bandages securely in place without adhesive tape. The breathable material provides a gentle, yet firm hold, reducing the risk of dressing displacement while ensuring maximum comfort.

KEY FEATURES

- No adhesive – protects sensitive skin
- Pre-cut and ready-to-use for fast application
- High elasticity for a secure and comfortable fit
- Reduces dressing displacement and wound exposure



CAREFIX® HEADBAND

SECURE FIXATION FOR HEAD, EAR & NECK DRESSINGS

Carefix® Headband is an elastic and comfortable fixation solution designed to secure surgical dressings and bandages after head, ear, or neck surgeries. It eliminates the need for adhesives, reducing the risk of skin irritation and ensuring a stable fit even for prolonged use.

KEY FEATURES

- Elastic and breathable fabric for enhanced comfort
- Easy and quick application without additional tools
- Pre-cut and ready-to-use to reduce waste
- Helps protect surgical sites and dressings from displacement

APPLICATIONS

- Fixation of post-operative dressings on the head, ears, and neck
- Protection and stabilization of wound care bandages
- Suitable for post-surgical and trauma care



CAREFIX® HEAD FIXATION

TAPE-FREE FIXATION FOR HEAD DRESSINGS

Carefix® Head Fixation is designed for patients requiring secure head dressing fixation without adhesives. Its elastic and soft material conforms to the head shape, ensuring a stable and comfortable fit while keeping surgical bandages and wound dressings securely in place.

KEY FEATURES

- No need for adhesive tapes, reducing skin irritation
- High elasticity ensures a snug and secure fit
- Pre-cut and ready-to-use for fast and efficient application
- Helps maintain wound protection, reducing the risk of infection

APPLICATIONS

- Fixation of surgical dressings on the head without tape
- Protection of post-operative wounds and sensitive areas
- Ideal for both clinical and at-home wound care



CAREFIX® FIXATION SOLUTIONS

ADDITIONAL PRODUCT INFORMATION

Composition & washing instructions

79% Polyamide, 21% Elastane



Sizes available

S-L

Packaging

Dispenser boxes & plastic bag

Intended use

Elastic dressing fixation solutions used to support and secure medical dressings or devices on various body parts, ensuring optimal wound care and healing.



WE'RE CERTIFIED

TYTEX IS CERTIFIED ACCORDING TO

